HUMANTECHNIK

lisa

D

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Seite 2

combi | *Sender*

GB

Operating Instructions

Page 8

combi | *transmitter*

F

Mode d'emploi

Page 12

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NL

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Pagina 18

combi | *Zender*

I

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Pagina 22

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

E

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Página 28

combi | *Emisor*

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Bestell-Nr.:

A-2950-0 (Länge 5 Meter)

A-2951-0 (Länge 10 Meter)

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

A-2951-0 (length 10 m)

Plug to connection sockets
combi I

Congratulations on purchasing the *lisa Combi I*. We hope it will give you many years of reassuring operation. Please read through these operating instructions carefully. They tell you how to set up the system correctly and explain all the system's features to you.

***lisa* system introduction**

An alerting system requires at least one transmitter and a receiver. The transmitter captures sound signals and transforms them into digital signals that are passed to a receiver by the existing mains electricity supply. The receiver then transforms these digital signals into light signals.


Channel selection

A digital transmission system ensures the correct recognition of each of the digital signals from each transmitter. The system has 10 different selectable channels, 0-9. The factory setting for all *lisa* devices is channel 0. Only where more than one installation is required in the same building will it be necessary to select a different channel for each system.

To select a channel, turn the rotary switch on the underside of the transmitter to the required channel number using a small screwdriver.

Example:

Set channel number 0

Personal Call button

Pressing this button causes a certain flash sequence on each receiver. It is possible to call people in other rooms this way. This button also allows the *lisa* system to be tested without any connecting cables.

Duration of flashing signals

The factory setting for the duration of the flash is 28 seconds for all *lisa* receivers. However, you can reduce this to 8 seconds if you wish.

To change the flash duration for a particular transmitter, remove the transmitter from the mains socket and, while pressing in the Personal Call button, plug the transmitter back into the socket. To switch back to 28 seconds remove the transmitter from the socket and plug back in again this time without pressing the Personal Call button.

Doorbell connection

Plug the transmitter into a convenient mains socket. Plug the 2.5 mm jack plug of the doorbell cable into the socket at the underside of the transmitter. Lay the cable to the bell push or door chime. Connect both wires to the screw contacts of the bell push or door chime.

Warning

The Transmitter is only designed to work with low voltage doorbells.

Do not connect the thin doorbell cable to the mains supply. This could be extremely dangerous for the user and may destroy the transmitter.

The doorbell voltage must be at least 3 Volts and must not exceed 15 Volts. Both AC and DC can be used.

Pressing your doorbell will now cause each receiver to flash with the corresponding number of flashes.

A second doorbell may be connected into the remaining doorbell socket. The receiver will flash differently for doorbell 1 and doorbell 2.

Installation to an intercom system

If your intercom is capable of differentiating between 2 different door chimes (e.g. the door to your flat and the entrance door to the building), you can convey the 2 door chimes to the transmitter using 2 connection cables. The correct assignment of the connecting cables to the intercom system is very important. It is recommended that you contact a qualified electrician and take note of the previous paragraph on connecting doorbells.

Failure causes

- You did not press the doorbell for long enough. Ring the bell for longer.
- If your doorbell system uses direct current voltage, you have not connected the wires with the correct polarity. Swap over the wires.

- The doorbell voltage of your system is too low. Either contact an electrician to have the doorbell system converted or use the *Doorbell Transmitter - Acoustic*.

Telephone connection

Plug the end of the cable with the transparent clip-in connector (6-pin plug) into the respective socket at the back of the transmitter. Plug the other end into the telephone socket.

If the required telephone socket is in use plug a 3-way telephone adaptor into the socket and plug your telephone into one of the sockets on the adaptor.

Warning:

Only use the supplied telephone cable otherwise correct operation can not be guaranteed.

The *lisa Combi I* is only designed for connection to analogue systems. If a connection is required to a digital system (e.g. ISDN or part of a PBX exchange system) please use a Telephone Transmitter – Acoustic or the *Universal lisa*.

Important note

All *lisa* devices connected to the mains supply have a low power consumption which may

lead to the device feeling warm. This is normal and does not interfere with the normal function of the devices.

These signalling systems have a high operational safety. Should any faults occur in spite of proper installation and use, please contact your dealer.

High and low voltage spikes (e.g. during a thunderstorm) can occur on the mains supply.

The *lisa* devices have an integrated electronic controller which may stop working after these conditions. If this happens simply unplug each of the devices and plug them in again.

As a precaution remove all devices from their sockets during a thunderstorm.

The use of halogen lamps or other equipment (e.g. Microwave Ovens) may, in rare cases, interfere with the *lisa* signals. In these cases the use of mains filters or suppressors is advisable. Please ask your dealer.

Maintenance and care

All *lisa* devices are maintenance free. To clean, remove from their mains sockets and wipe over with a moistened cloth. Never use alcohol, thinner or other organic materials.

lisa devices should not be exposed to excessive direct sunlight, heat, thick dust, humidity or mechanical shock.

Warranty

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.



The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people who are not authorised to do so (destruction of the seal on the unit). Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer.

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Humantechnik GHG AG