*lisa*

**D**

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Seite 2

combi III · Sender

**GB**

**Operating Instructions**

Page 8

combi III · transmitter

**F**

**Mode d'emploi**

Page 12

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**NL**

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Pagina 18

combi III · zender

**I**

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Pagina 22




combi III · trasmettitore

**E**

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Página 28

combi III · emisor

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





Plug to connection to combi III

Microphone to stick over the acoustic output

Congratulations on purchasing the *lisa Combi III*. We hope it will give you many years of reassuring operation. Please read through these operating instructions carefully. They tell you how to set up the system correctly and explain all the system's features to you.

### **lisa system introduction**

An alerting system requires at least one transmitter and a receiver. The transmitter captures sound signals and transforms them into digital signals that are passed to a receiver by the existing mains electricity supply. The receiver then transforms these digital signals into light signals.


### **Channel selection**

A digital transmission system ensures the correct recognition of each of the digital signals from each transmitter. The system has 10 different selectable channels, 0-9. The factory setting for all lisa devices is channel 0. Only where more than one installation is required in the same building will it be necessary to select a different channel for each system.

To select a channel, turn the rotary switch on the underside of the transmitter to the required channel number using a small screwdriver.

### **Example:**

Set channel number 0

### **Personal Call button**

Pressing this button causes a certain flash sequence on each receiver. It is possible to call people in other rooms this way. This button also allows the lisa system to be tested without any connecting cables.

### **Duration of flashing signals**

The factory setting for the duration of the flash is 28 seconds for all lisa receivers. However, you can reduce this to 8 seconds if you wish.

To change the flash duration for a particular transmitter, remove the transmitter from the mains socket and, while pressing in the Personal Call button, plug the transmitter back into the socket. To switch back to 28 seconds remove the transmitter from the socket and plug back in again this time without pressing the Personal Call button.

### **Doorbell connection**

Plug the transmitter into a convenient mains socket. Plug the doorbell microphone lead into the respective socket at the back of the

transmitter. Lay the cable to the door chime and attach the microphone to the door chime by peeling off the paper circle to reveal the sticky ring beneath.

For optimal operation of the transmitter the sensitivity may be adjusted so that only the doorbell activates the transmitter and not background noise.

### Setting the sensitivity

1. Set the sensitivity controller to its lowest level (turn to the anticlockwise to stop)
2. Activate the doorbell and gradually increase the sensitivity of the transmitter (by turning the dial clockwise) until the transmitter reacts to the doorbell.
3. If your intercom system operates with different signal tones, make sure that both signal tones are reliably detected.

### Telephone connection

Plug the microphone cable into the socket marked at the underside of the transmitter. Lay the cable to the telephone. Peel off the paper circle on the microphone to reveal the sticky ring beneath. Attach the microphone to the area of the telephone where the ringer is situated. Adjust the volume of the ringer so that the transmitter reacts to the ringing.

Note: If the telephone handset is

replaced too harshly it may activate the telephone transmitter.

### Important note

All lisa devices connected to the mains supply have a low power consumption which may lead to the device feeling warm. This is normal and does not interfere with the normal function of the devices.

These signalling systems have a high operational safety. Should any faults occur in spite of proper installation and use, please contact your dealer.

High and low voltage spikes (e.g. during a thunderstorm) can occur on the mains supply. The lisa devices have an integrated electronic controller which may stop working after these conditions. If this happens simply unplug each of the devices and plug them in again.

**As a precaution remove all devices from their sockets during a thunderstorm.**

The use of halogen lamps or other equipment (e.g. Microwave Ovens) may, in rare cases, interfere with the lisa signals. In these cases the use of mains filters or suppressors is advisable. Please ask your dealer.

### Maintenance and care

All lisa devices are maintenance free. To clean, remove from their mains sockets and wipe over with a

moistened cloth. Never use alcohol, thinner or other organic materials. lisa devices should not be exposed to excessive direct sunlight, heat, thick dust, humidity or mechanical shock.

**Note:** This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

#### **Environmental note/Waste disposal**

Disposal of used electrical and electronic equipment (applicable in the countries of the European Union and in other European countries which operate a separate collection system for these devices).

The symbol on the product or its packaging indicates that this product should not be treated as normal household waste. It should be handed in to a collection center for recycling electrical and electronic equipment.

By playing your part in ensuring the correct disposal of this product, you will be safeguarding the environment and the health of other people. The environment and health are jeopardized by incorrect disposal. Material recycling helps to reduce the consumption of raw materials. You can obtain further information

about recycling this product from your local council, from local waste disposal operators and from the shop where you bought the product.

#### **Warranty**

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people who are not authorised to do so (destruction of the seal on the unit).




Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer.

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### **Humantechnik GHL AG**