HUMANTECHNIK

*lisa*

**DE**

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Seite 2

Funk-Tischblinklampe

**EN**

**Operation Instructions**

Page 8

RF table blink lamp

**FR**

**Mode d'emploi**

Page 14

Lampe-flash de table radio *lisa*

**NL**

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Pagina 20

RF tafelknipperlamp

**IT**

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Pagina 26

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



## Side view right

## Side view left

Congratulations on purchasing your »RF table blink lamp«. You have chosen a modern and reliable receiver for our lisa RF system. Please read this operating instruction carefully to be able to set up the unit correctly and to familiarise yourself with all of the systems features.

### **Standard components**

Please check if all of the following components are included:

- »RF table blink lamp«
- 12 V plug-in PSU
- 9 V rechargeable battery
- Operating Instructions
- Warranty card

Should any parts be missing, please immediately contact your dealer or the manufacturer directly.

### **Function principle**

The »RF table blink lamp« picks up signals from different radio transmitters and converts them into light signals.

The »RF table blink lamp« has a very powerful light-emitting diode (LED) which brings a signal to the user's attention. Further additional LEDs indicate the nature of the incoming signals.

### **Controls**

#### **Battery/Depleted battery**

If the »Battery/Dead battery« indicator lights up constantly, this means that the battery or rechargeable battery which is inserted is almost fully discharged. A regular battery must be replaced or the rechargeable battery recharged. The »Depleted battery« indicator might flash briefly in the flashing mode; this does not indicate that the battery is low.

### **Indicator LEDs**

The 6 coloured indicator light-emitting diodes (LED) help you to differentiate between the individual transmitter signals. They light up for a total of 40 seconds.

### **lisa-Channel selector**

You can select the radio channel using the channel selector. Make sure all devices in your system are set to the same channel.

### **ON-OFF switch**

Set the on/off switch to »ON«. The »RF table blink lamp« is now ready to operate.

After the unit has been switched on, each of the coloured indicator LEDs will light up briefly in succession. This forms part of a self-test of the unit and lets you know that it is functioning correctly.

### **OK button**

Press the OK button to acknowledge incoming signals. This will switch off the bright blinking LED. The coloured indicator LED lights will switch off automatically after a total of 40 seconds.

Furthermore, you can press the OK button to call up the most recent signal again (the function works for up to 3 minutes after reception of the signal). The strong white blinking LED signal will not be repeated.

### **Note (special accessories)**

Vibrating pillow or other modules:

- Acoustic module MA-1
- Flash module MF-1
- Switch module MS-1

can be connected to the corresponding jacks of the »RF table blink lamp«.

## Setting up the system

### 1a. Mains operation

Insert the plug at the end of the 12 V power supply into the corresponding jack on the rear side of the »RF table blink lamp«. Then plug the power supply unit into the normal 230 V socket in your house.

### Power failure backup

Please note that the device is designed for mains operation. The rechargeable battery or normal battery are used as backup only, in the event of a power failure.

### 1a. Inserting the rechargeable battery (included)

There is a battery compartment on the back of the »RF table blink lamp«. Open it by sliding the cover of the battery compartment to the right of the receiver using your thumb. Connect the battery plug with the outlet, insert the rechargeable battery into the compartment and close the device. When the battery is connected, the receiving signals may be triggered once.

### 1c. 1b. Operation with a 9 V battery (not included)

The »RF table blink lamp« can also be operated using a regular 9 V block battery, which is widely available. Open the battery compartment as described in 1a. Connect the battery to the connector plug and insert the battery into the battery compartment. Close the device. is connected to the mains power supply.


**Prior to using the rechargeable battery for the first time, it must be fully charged for 12 hours.**

**For rechargeable battery operation, only use and recharge the special rechargeable battery A-2995-0. Rechargeable batteries cannot be recharged using the 9 V block battery connection!**

## 2. Setting the correct radio channel

A special digital transmission process ensures that the radio signals are picked up without errors. This process also makes it possible to set 10 different radio channels. When supplied from the factory, all »RF table blink lamp« are set to channel 0. It is only necessary to set units to other channels if you have several different users operating RF lisa systems in the same building. Set the required channel by using a small screwdriver to turn the channel selector switch, which is located on the right side of the unit, to the required channel number. All transmitters and receivers must be set to the same channel.









Example: Set channel number 0

## Transmission signal detection

The »RF table blink lamp« converts the signals sent by the *lisa* RF transmitter into different sequences of blinking and lamp signals.

At the same time, it is also possible to see which *lisa* signal has been sent because the corresponding coloured LED will either flash or light up for 40 seconds.

**P12** Pause 1,2 seconds

**P14** Pause 1,4 seconds

**P20** Pause 2,0 seconds

**B** Blink

**P30** Pause 3,0 seconds

**P35** Pause 3,5 seconds

**P40** Pause 4,0 seconds

**W** Repeat

The 6 coloured indicator LEDs remain on for 40 seconds.

## Radio range

Note the following points when using the »RF table blink lamp«:

The radio signals are designed to have a range of 80 meters under perfect conditions. The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation.
- Interference from TVs, computer monitors and mobile phones. Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavourable distance from surfaces which reflect radio waves, for example floors and walls. This may lead to the signal wave being weakened or even cancelled out by the reflected wave.
- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

## Maintenance and care

The »RF table blink lamp« does not require any maintenance. If the unit does become dirty, simply wipe it clean with a soft, damp cloth. **Note: Disconnect the »RF table blink lamp« from the power supply!**

Never use spirits, thinners or other organic solvents. Do not set up the »RF table blink lamp« where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

**Note:** This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Please make sure that the batteries are not exposed to excessive heat, such as sunlight, fire or anything similar.

## Warranty

The »RF table blink lamp« is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge. It is recommended that you send in the product in its original packaging, so keep the packaging for the duration of the warranty period.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by unauthorised persons (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

**Always specify the product number in any event.**

**Disposal** of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

### Technical data

Power supply:	12 V-power supply unit at 230 V, 50 Hz via a 9 V rechargeable battery A-2995-0 (included) or 9 V block battery (not included)	
Battery standby time:	Approx. 1-2 days (at 5 events a day)	
Rechargeable battery standby time:	Approx. 10 hours (at 5 events a day)	
Frequency:	868.35 MHz	
Operating range:	0° C to 40° C	
Blink signaling:	1 LED high power with 2.5 W of power	
Battery/Dead battery indicator:	1 LED orange	
Visual signals:	1 LED red (fire)	1 LED orange (telephone)
	1 LED blue (water)	1 LED yellow (pager)
	1 LED white (baby phone)	1 LED green (doorbell)
	all 6 LEDs (alarm)	
Height:	65 mm	
Width:	93 mm	
Depth:	150 mm	
Weight:	220 g	

This device satisfies the following EC directives:


- 2011 / 65 / EC RoHS directive
- 2012 / 19 / EC WEEE directive
- 2014 / 53 / EC RED directive


Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at [www.humanttechnik.com](http://www.humanttechnik.com).



Technical specifications subject to change without notice.

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**Sarabec Ltd**