lisa

D

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Seite 2

Funk-Blinklampe

GB

Operating Instructions

Page 10

RF blink lamp

F

Mode d'emploi

Page 18

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NL

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Pagina 26

Radio-knipperlamp

I

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Pagina 34


Ricevitore radio luminoso

E

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Página 42

Radorreceptor de luz intermitente

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
Congratulations on purchasing your *RF blink lamp*. We hope it will give you many years of reassuring operation. You have purchased a modern and reliable system.

Please read through these operating instructions carefully. They tell you how to set up the system correctly and explain all the system's features to you.

Standard components

Please check if all following components are included:

- *RF blink lamp*
- 9 V block battery
- 9 V plug-in PSU
- Stand
- Operating Instructions
- Warranty card

If any parts are missing please immediately contact your dealer or the manufacturer directly.


Functional principle

The *RF blink lamp* receives signals from a variety of radio transmitters and converts them into light signals. The *RF blink lamp* has a very high-power light-emitting diode (LED), which brings a signal to the user's attention, and additional LEDs which indicate the nature of the incoming signals.

Setting up the system

1a. Inserting the battery

There is a battery compartment on the back of the *RF blink lamp*. Open it by sliding the cover of the battery compartment to the right of the receiver using your thumb. Connect a 9 V block battery to the connector plug, insert the battery into the battery compartment and close the device. When the battery is connected, the receiving signals may be triggered once.

1b. Operation with a rechargeable battery (A-2995-0)

The *RF blink lamp* can also be operated using a rechargeable battery (A-2995-0), which is available as a special accessory. Open it by sliding the cover of the battery compartment to the right of the receiver using your thumb. To connect the rechargeable battery, insert the battery plug into the outlet in the battery compartment.

The separate rechargeable battery connection makes sure that a 9 V block battery will not be recharged by mistake. The rechargeable battery is recharged when the *RF blink lamp* is connected to the mains power supply.

Prior to using the rechargeable battery for the first time, it must be recharged for 12 hours.

For rechargeable battery operation, use and recharge only the special rechargeable battery A-2995-0.

Rechargeable batteries cannot be recharged using the 9 V block battery connection!

2. Setting the correct radio channel

A special digital transmission process ensures that the radio signals are picked up without errors. This process also makes it possible to set 10 different radio channels. When supplied from the factory, all *RF blink lamps* are set to channel 0.

It is only necessary to set units to other channels if you have several different users operating units in the same building. Set the required channel by using a small screwdriver to turn the channel selector switch, which is located towards the top of the battery compartment, to the required channel number.

All transmitters and receivers must be set to the same channel.

Fig. setting the correct radio channel.

Operation

1. Switching on the unit

Set the on/off switch to on. The *RF blink lamp* is now ready to operate. After the unit has been switched on, each of the indicator lamps lights up briefly in succession and the high-power light-emitting diode blinks several times. This forms part of a self-test of the unit and lets you know that it is functioning correctly.


2. Transmission signal detection

The *RF blink lamp* converts the various transmission signals into different sequences of blinking and lamp signals.

Telephone: Telephone LED lights

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P14', then another box labeled 'B', followed by a box labeled 'P30', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.

Doorbell I: Doorbell LED lights

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P14', then another box labeled 'B', followed by a box labeled 'P14', then another box labeled 'B', followed by a box labeled 'P14', then another box labeled 'B', followed by a box labeled 'P40', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.

Doorbell II: Doorbell LED flashes

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P40', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.

Babymonitor: Babymonitor LED lights

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P35', then another box labeled 'B', followed by a box labeled 'P40', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.


Person call: Person call LED lights

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P12', then another box labeled 'B', followed by a box labeled 'P12', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.


Alarm: All LEDs flash

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P14', then another box labeled 'B', followed by a box labeled 'P14', then another box labeled 'B', followed by a box labeled 'P14', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.

Water alarm: Water alarm LED lights

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P12', then another box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P12', then another box labeled 'B', followed by a box labeled 'P12', then another box labeled 'B', followed by a box labeled 'P12', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.

Fire alarm: Fire alarm LED lights

Blink rhythm  A horizontal timeline showing a sequence of events: a box labeled 'B', followed by a box labeled 'P12', then another box labeled 'B', followed by a box labeled 'P20', then another box labeled 'B', followed by a box labeled 'P12', then another box labeled 'B', followed by a box labeled 'P20', and finally a box labeled 'W'. Each box is connected to the next by a horizontal line.

P12 pause 1,2 seconds

P14 pause 1,4 seconds

P20 pause 2,0 seconds

B flash

P30 pause 3,0 seconds

P35 pause 3,5 seconds

P40 pause 4,0 seconds

W repeat

The indicator LEDs remain on for 40 seconds.

Battery/rechargeable battery empty

If the »battery/rechargeable battery empty-display« lights up constantly, this means that the battery or rechargeable battery used is almost fully discharged. A regular battery must be replaced or a rechargeable battery recharged.

The »battery/rechargeable battery empty-display« might flash briefly in the flashing mode; this does not indicate that the battery is low.

Acoustic alarm

The acoustic alarm can be switched on or off using the on/off button.

Reset button

Press the reset button to acknowledge incoming signals. This will switch off the flashing light and the acoustic signal.

Retrieving the last signal

The *RF blink lamp* stores any signal it has received for approx. 3 minutes. During this time, you can repeat the signal by pressing the reset button again. The signal will then be indicated by the signal lights. The flashing light and the acoustic signal will not be repeated.

Mains operation

You can also operate your *RF blink lamp* via the plug-in power supply unit. To connect the PSU, insert the connector on the end of the PSU cable into the corresponding socket on the left-hand side of the device.

Then plug the power supply unit into the normal 230 V socket in your house.

Recharging the A-2995-0 rechargeable battery

If there is a rechargeable battery in the *RF blink lamp* while the lamp is plugged into a power supply, the battery will be recharged.

The recharging cycle begins automatically. After 12 hours, the device will switch to conservation charging. The device can remain plugged into the regular power supply. The *lisa* device's reception functions will continue to operate while recharging is in progress.

Radio range

The radio signals are designed to have a range of 80 meters under optimum conditions.

The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation
- Interference from TVs, computer monitors and mobile phones.
Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavourable distance from surfaces which reflect radio waves, for example floors and walls.
This may lead to the signal wave being weakened or even cancelled out by the reflected wave.
- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

Maintenance and care

The *RF blink lamp* does not require any maintenance. If the unit does become dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents.

Do not set up the *RF blink lamp* where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Please make sure that the batteries are not exposed to excessive heat from such as sunlight, fire or anything similar.

Warranty

The *RF blink lamp* is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly. This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away. The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system). The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.


Technical data

Power supply:	9 V block battery, via a plug-in PSU or rechargeable battery A-2995-0 (special accessory)
Battery life:	Approx. 1-2 weeks (assuming 5 events per day)
Frequency:	868.35 MHz
Operating range:	0° C – 40° C
Flash signaling:	1 LED high power with 2.5 W of power
Battery-low indicator:	1 LED orange
Visual signals:	2 LED red (fire/water alarm) 2 LED yellow (telephone/doorbell signal) 2 LED orange (baby phone/pager) All 6 LEDs (alarm)
Height (including the antenna):	147 mm (201 mm)
Width:	50 mm
Depth:	41 mm
Weight (including battery):	150 g

This product complies with European Union norms.
You can obtain EC compliance declarations from your dealer, or directly from the manufacturer of this product.

Technical specifications subject to change without notice.

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Humantechnik GHL AG