HUMANTECHNIK

lisa

DE

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Seite 2

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EN

Operation Instructions

Page 8

Acoustic RF receiver

FR

Mode d'emploi

Page 14

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NL

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Pagina 20

Akoestische radio-ontvanger

IT

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Pagina 26

Radoricevitore acustico


ES

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
Página 32

Radorreceptor acústico

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





Side view right

Side view left

Congratulations on purchasing your *acoustic RF receiver*. We hope it will provide you with many years of service. You have chosen a modern and reliable system.

Please read through these operating instructions carefully to be able to start the unit correctly and to familiarise yourself with all of the systems features.

Standard components

Please check if all following components are included:

- *Acoustic RF receiver*
- 9 V plug-in PSU
- Operating Instructions
- Warranty card

If any parts are missing please immediately contact your dealer or the manufacturer directly.

Function

The *acoustic RF receiver* picks up signals from different lisa RF radio transmitters and converts them into audible alerts.

In addition, 6 event indicator LEDs help you to identify the various alerts.

Controls

Battery/Dead battery

If the battery status indicator lights up constantly, this means that the battery or rechargeable battery used is almost depleted. A regular battery must be replaced or a rechargeable battery recharged. The battery status indicator might flash briefly in the flashing mode; this does not indicate that the battery is low.

Event indicator LEDs

The 6 indicator lamps help you to differentiate between the individual transmitter signals. They light up for a total of 40 seconds.

***lisa*-Channel selector**

You can select the radio channel using the channel selector.

ON-/OFF switch

Set the on/off switch to »ON«. The *acoustic RF receiver* is now ready to operate. After the unit has been switched on, each of the indicator lamps lights up briefly in succession. This forms part of a self-test of the unit and lets you know that it is functioning correctly.

OK button

Press the OK button to acknowledge incoming signals. This will switch off the acoustic signal. The indicator lamps light and go out automatically after a total of 40 seconds. Furthermore, you can press the OK button to call up the most recent signal again (the function works for up to 3 minutes after receipt of the signal).

The acoustic signal will not be repeated.

Setting of the volume

Volume: continuously adjustable
(max. 100 dB)

Note (special accessory)

Vibrating pillow or other modules:

- Flash module MF-1
 - Switch module MS-1
- can be connected to the corresponding jacks of the *acoustic RF receiver*.

The following audio devices can be connected to the 3.5 mm audio out jack:

- Infrared transmission systems
- Radio transmission systems
- Induction loop transmission systems
- Hi-fi systems or active speaker systems

Setting up the system

1a. Mains operation

Insert the plug at the end of the 9 V power supply into the corresponding jack on the rear side of the *acoustic RF receiver*. Then plug the power supply unit into the normal 230 V socket in your house.

Power failure backup

Please note that the device is designed for mains operation. The battery and rechargeable battery are used as backup only in the event of a power failure.

1b. Inserting the battery

The *acoustic RF receiver* can also be powered by a 9 V block battery (not included with the device). There is a battery compartment located on the rear side of the *acoustic RF receiver*. Open it by sliding the cover of the battery compartment to the right of the receiver using your thumb. Connect a 9 V block battery to the connector plug, place the battery into the recess and close the receiver again. When the battery is connected, the receiving signals may be triggered once.

1c. Inserting the rechargeable battery

The *acoustic RF receiver* can also be powered by a rechargeable battery (A-2995-0), which is available as a special accessory. Open it by sliding the cover of the battery compartment to the right of the receiver using your thumb. To connect the rechargeable battery, insert the battery plug into the outlet in the battery compartment. The separate rechargeable battery connection makes sure that a 9 V block battery will not be recharged by mistake. The rechargeable battery is

recharged when the *acoustic RF receiver* is connected to the mains power supply.

Prior to using the rechargeable battery for the first time, it must be recharged for 12 hours.


For rechargeable battery operation, use and recharge only the special rechargeable battery A-2995-0. Rechargeable batteries cannot be recharged using the 9 V block battery connection!

2. Setting the correct radio channel

A special digital transmission process ensures that the radio signals are picked up without errors. This process also makes it possible to set 10 different radio channels. When supplied from the factory, all *acoustic RF receivers* are set to channel 0. It is only necessary to set units to other channels if you have several different users operating units in the same building. Set the required channel by using a small screwdriver to turn the channel selector switch, which is located on the right side of the unit, to the required channel number.

All transmitters and receivers must be set to the same channel.

Example: Set channel number 0

Transmission signal detection

The *acoustic RF receiver* converts the signals sent by the *lisa* RF transmitter into the different sounds (telephone ringing, doorbell, baby cries, etc.).

At the same time, it is also possible to see which *lisa* signal has been sent because the corresponding LED will either flash or light up for 40 seconds.

Telephone: Telephone LED lights up in orange and you can hear the sound of a telephone ringing.

Door I: Door LED lights up in green and you can hear the sound of a doorbell.

Door II: Door LED flashes in green and you can hear the sound of a different doorbell to the first one.

Babycall: Babymonitor LED lights up in white and you can hear the sound of a baby crying.

Person call: Person call LED lights up in yellow and you can hear the pager signal.

Alarm: All 6 LEDs flash and you can hear an alarm signal.

Water alarm: Water-alarm LED lights up in blue and you can hear an alarm signal.

Fire alarm: Fire-alarm LED lights up in red and you can hear an alarm signal.

Transmission range

Note the following points when using the *acoustic RF receiver*:

The radio signals are designed to have a range of 80 meters under optimum conditions. The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation.
- Interference from TVs, computer monitors and mobile phones. Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavourable distance from surfaces which reflect radio waves, for example floors and walls. This may lead to the signal wave being weakened or even cancelled out by the reflected wave.
- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

Maintenance and care

The *acoustic RF receiver* does not require any maintenance. If the unit does become dirty, simply wipe it clean with a soft, damp cloth.

Note: remove the *acoustic RF receiver* from mains socket!

Never use spirits, thinners or other organic solvents. Do not set up the *acoustic RF receiver* where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Please make sure that the batteries are not exposed to excessive heat from such as sunlight, fire or anything similar.

Warranty

The *acoustic RF receiver* is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge. It is recommended that you send in the product in its original packaging, so keep the packaging for the duration of the warranty period.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

Technical data

Power supply:	9 V-power supply unit at 230 V, 50 Hz via a 9 V block battery or rechargeable battery A-2995-0 (special accessory)
Battery standby time: Rechargeable battery standby time:	Approx. 2-3 days (at 5 events a day) Approx. 10 hours (at 5 events a day)
Frequency:	868.35 MHz
Temperature range:	0° C to 40° C
Battery/Dead battery indicator:	1 LED orange
Visual alerts:	1 LED red (fire alert) 1 LED blue (water alert) 1 LED orange (telephone alert) 1 LED green (doorbell alert) 1 LED white (baby call) 1 LED Gelb (person call) all 6 LEDs active: Alarm
Height:	65 mm
Width:	93 mm
Depth:	134 mm
Weight:	220 g

This device complies with the following EC directives:

2011 / 65 / EC + 2015 / 863 / EC RoHS-directive


2012 / 19 / EC WEEE-directive

2014 / 53 / EC RED-directive


Compliance with the directives listed above is confirmed by the CE seal on the device.
EC compliance declarations are available on the Internet at www.humanteknik.com.



Technical specifications subject to change without prior notice.

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Sarabec Ltd