HUMANTECHNIK

lisa
TV-Link

D

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Seite 2

lisa RF TV-Link

GB

Operating Instructions

Page 8

lisa RF TV-Link

F

Mode d'emploi

Page 14

lisa RF TV-Link

NL

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Pagina 20

lisa RF TV-Link

I

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Pagina 26

lisa RF TV-Link


E

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Página 32


lisa RF TV-Link

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Rückansicht

lisa-
Kanalwahl-
schalter

Audio
in

Mono
Stereo

Audio
out

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
Front view

Back view

Congratulations on purchasing your *lisa* RF TV-Link. We hope it will provide you with many years of service. You have chosen a modern and reliable system.

Please read the operating instructions through carefully to be able to start the unit correctly and to familiarise yourself with all of the systems features.

Standard package

Please make sure that all the following parts have been included in the package:

- *lisa* RF TV-Link
- 230 V mains power cord
- TV connection cord set
- Operating instructions
- Guarantee card

If any parts are missing please immediately contact your dealer or the manufacturer directly.

How it works

The *lisa* RF TV-Link picks up signals from different *lisa* transmitters, converts them into the corresponding sound, which is then transmitted to the TV hearing-aid system connected, e.g. from the *InfraLight* series. 6 signal lamps (LEDs) also visually indicate the different transmission signals.

Display and operating elements/signal lamps

Signal lamps

The 6 signal lamps (LEDs) visually indicate which *lisa* signal has been sent. The display switches off after 40 seconds.

***lisa*-channel selector**

Use this switch to select the *lisa* RF TV-Link radio channel.

Mono/stereo switch

If the *lisa* RF TV-Link is operated on a mono sound source, this switch must be set to Mono.


Installation

1. Setting the correct radio channel

A special digital transmission process ensures that the radio signals are picked up without errors. This process also makes it possible to set 10 different radio channels. When supplied from the factory, all *lisa* units are set to channel 0. It is only necessary to set units to other radio channels if there are several systems being operated at the same time in the same building.


Set the required channel by using a small screwdriver to turn the channel selector switch, which is located towards the top of the battery compartment, to the required channel number.

All transmitters and receivers must be set to the same channel.

Example: Set channel number 0

2. Installation of *lisa* RF TV-Link

Important Note

The *lisa* TV link is compatible with TV listening systems in the *InfraLight* series or

comparable infrared TV listening systems and popular loop amplifiers (e.g. LA-215).

When using a RadioLight or comparable wireless TV listening system operating on a frequency of 863 HMz, you will achieve the best results by **keeping a minimum distance of 2 meters from the *lisa* TV link**. This is because the signals from the wireless TV listening system and the *lisa* may interfere with one another.

Operation

1. Signal sequence

Once a *lisa* radio signal is detected by the *lisa* RF TV-Link, this device converts the signal into the corresponding sound and transmits it to the TV hearing-aid system connected. The *lisa* RF TV-Link transmits the sound for approx. 30 seconds, during this time, the television volume is turned down.

2. Transmission signal detection

The *lisa* RF TV-Link converts the signals sent by the *lisa* RF transmitter into the different sounds (telephone ringing, doorbell, baby cries, etc.).

At the same time, it is also possible to see which *lisa* signal has been sent because the corresponding LED will either flash or light up for 40 seconds.

Fire alarm: Fire-alarm LED lights up and you can hear the sound of a fire-engine siren.

Water alarm: Water-alarm LED lights up and you can hear the sound of water running.

Telephone: Telephone LED lights up and you can hear the sound of a telephone ringing in your TV hearing-aid system.

Door I: Door LED lights up and you can hear the sound of a doorbell.

Door II: Door LED flashes and you can hear the sound of a different doorbell to the first one.

Babyphone: Baby LED lights up and you can hear the sound of a baby crying.

Pager: Pager LED lights up and you can hear the pager signal.

Alarm: All 6 LEDs flash and you can hear an alarm signal.

Radio range

Note the following points when using the *lisa* RF TV-Link:

The radio signals are designed to have a range of 80 meters under optimum conditions.

The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation
- Interference from TVs, computer monitors and mobile phones. Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavourable distance from surfaces which reflect radio waves, for example floors and walls. This may lead to the signal wave being weakened or even cancelled out by the reflected wave.
- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

Maintenance and care

The *lisa* RF TV-Link does not require any maintenance. If the unit does become dirty, simply wipe it clean with a soft, damp cloth. **Note:** **remove the *lisa* TV link from mains socket!** Never use spirits, thinners or other organic solvents. Do not set up the *lisa* RF TV-Link

where it will be exposed to full sunlight for long periods.

In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is not protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Warranty

The *lisa* RF TV-Link is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system). The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.


Specifications

Power supply:	230 V, 50 Hz
Frequency:	868.35 MHz
Visual signals:	2 LEDs, red (fire/water alarm) 2 LEDs, yellow (telephone/doorbell signal) 2 LEDs, orange (babyphone/pager)
Height:	41 mm
Width:	117 mm
Depth:	133 mm
Weight:	Approx. 300 g
Temperature range:	0° C to 40° C


Technical specifications subject to change without notice.

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Humantechnik GHL AG