HUMANTECHNIK

signolux

DE

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Seite 2

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




③

Charging contacts On/Off Function button

Prop

Belt clip

Charging socket AUX-connector Delete button

Indicator
LEDs,
event
icons

Low
battery
indicator

Battery
compartment

signolux

Congratulations on purchasing your »signolux pager«!

Please read this manual carefully. It will tell you how to set up the system correctly and explain all the features it has to offer.

Contents of package

Please check if the following parts are included:

- »signolux-pager«
- neck cord
- 2 x 1.5 V AAA batteries
- User manual
- Warranty card

If any parts are missing, please contact your dealer or the manufacturer immediately.

Function

The portable »signolux-pager« is a receiver for the alerts of the various »signolux« transmitters and indicates these as vibrations.

6 indicator LEDs with symbols help to distinguish the alert source.

Controls and displays

Event indicator LEDs:

The 6 indicator lamps help you to differentiate between the individual transmitter signals. They light up for a total of 40 seconds.

Battery status indicator:

The receiver informs about the state of the batteries of the receiver as well as the batteries of the transmitters:

- A sufficient battery level of the receiver is indicated by the indicator LED flashing green every 5 to 6 seconds.
- A low battery level of the receiver is indicated by the indicator LED flashing red in short sequence (once per second).

Please replace the receiver batteries as soon as possible.

- A low battery level of a transmitter is indicated by the indicator LED shining together with the symbol of the according transmitter when it's signal is triggered.

Please replace the batteries of the according transmitter as soon as possible.

Function button:

When you feel the »signolux-pager« vibrating, you can stop the vibration by pressing the function button on the receiver.

Belt clip:

You can use the belt clip to attach the »signolux-pager« to your clothing. The belt clip can be exchanged with the neck cord (see image on page 17).

On/off switch:

Use this switch to turn the pager on and off.

Delete button:

This button serves to remove transmitters assigned to an indicator symbol (see section 4 „Deleting transmitters from receiver memory“).

USB-C charging socket:

There is a USB-C charging socket on the side of the receiver. If the receiver is equipped with rechargeable batteries (accessory: 2 x A-4966-0 rechargeable batteries, also included in the A-2910-0 charging station set), the receiver can also be charged using a standard USB charger with a USB-C cable instead of the charging station. During the charging process, the battery status indicator lights up red.

After the charging process has been completed, it lights up green.

The maintenance function of the charge control automatically keeps the batteries at full capacity. This means that the charging process can be active again at irregular intervals, with the battery status indicator lighting up red accordingly.

When operating with normal batteries this charging socket remains without function.

Setup

1. Inserting the battery

There is a battery compartment on the underside of the »signolux-pager«. Open it by sliding the cover of the battery compartment away from the receiver using both your thumbs.

Insert two 1.5 V AAA batteries and close the receiver again.

2. Switching the unit on

Set the on/off switch to »on«. The »signolux-pager« is now ready to operate. After the unit has been switched on, the battery status LED lights up for about a second.

Preset transmitter indication

Signalux devices can be used together immediately and without prior pairing when delivered.

Which event icon is displayed depends on the type of transmitter. This information can be found in the user manual of the according transmitter.

Optional manual pairing

Note 1:

If a transmitter is to be displayed with an event icon of your own choice, this can be done as described in the chapter Pairing transmitters and adapting them to the event icons.

However, this causes the alarm clock to lose the preset transmitter indication and all transmitters previously shown on the alarm clock must be individually paired.

Note 2:

Should there be repeated false alerts when using the preset transmitter indication, this may hint at a neighbouring Signalux installation causing interference. In this case, the problem can be solved by individually pairing all transmitter

and receiver units within the setup (see chapter 3 „*Pairing transmitters and adapting them to the event icons*“).

Reset to factory defaults

Delete the pairing of each event icon (see chapter 4 „*Deleting transmitters from the receiver memory*“).

When all previously assigned transmitters have been removed from the memory, the preset transmitter indication is restored and the transmitters are automatically displayed.

3. Pairing transmitters and adapting them to the event icons

A transmitter (e.g. doorbell pushbutton) can be assigned to each iconl. This assignment can be adapted by programming. Up to 8 transmitters can be assigned to each symbol. Please keep in mind that a transmitter can only be assigned if it is not yet known in the memory of the »signolux-pager«. If required, you might first have to delete the transmitter from the memory (see section 4.)

- a) Keep the function button pressed for some time. The symbols begin to light up in sequence. Let go of the function button as soon as the desired symbol is lit. (e. g. doorbell symbol)

- b) Trigger the signal of the transmitter you want to assign to the lit up symbol. (e.g. operate the doorbell button).

- c) Wait for about 20 seconds and test the function.

4. Deleting transmitters from the receiver memory

If you wish to remove the transmitters assigned to an icon from the receiver memory:

- a) Keep the function button pressed for some time. The icons begin to light up in sequence. Let go of the function button as soon as the desired icon is lit.
- b) Now press the function button and the delete button at the same time. All icons will briefly light up.

Please note that this will delete all transmitters assigned to the according symbol.

Icon description

doorbell (green)

telephone (orange)

person call (yellow)

babycry (white)

fire alarm (red)

alarm (red)

Checking the last alert, cancelling the current alert

By pressing the function button, the last registered alert can be displayed again.

You can also cancel an active alert by pressing this button.

Event identification

The various alert sources can be identified by the different vibration patterns which are unique to each icon:

Doorbell: Doorbell icon (green LED) is lit

Telephone: Telephone icon (orange LED) is lit

Person call: Person icon (yellow LED) is lit

Babycry: Baby icon (white LED) is lit

Fire alarm: Fire icon (red LED) is lit

Alarm: Alarm icon (red LED) is lit

Special accessories

The following modules (not included) can be plugged into to the AUX output connector:

- Vibrating pillow (A-3303-0)
- Flash module MF-1 (A-3320-0)
- Acoustic module MA-1 (A-3340-0)
- Vibrating contraption (A-3310-0)

Transmission range

The signal range of the »signolux« transmitters can be up to 200 metres under ideal conditions. The following conditions may result in a reduced range:

- Obstacles like building structures or vegetation.
- Interference from TVs, computer monitors and mobile phones. Under certain circumstances, these may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavourable distance from surfaces which reflect radio waves, for example floors and walls. This may lead to the signal wave being weakened or even cancelled out by the reflected wave.

- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Other devices operating at similar frequencies in the vicinity may also result in mutual interference.

Maintenance and care

All »signolux« devices are maintenance-free. If the units become dirty, simply wipe them clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents.

Please do not expose »signolux« devices to direct sunlight for a long time and please protect them from excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Please make sure that the batteries are not exposed to excessive heat from such as sunlight, fire or anything similar.

Warranty

All signalux devices are very reliable products. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly. This warranty covers the repair of the product and returning it to you free of charge.

It is recommended that you send in the product in its original packaging, so do not throw the packaging away while the warranty period lasts. The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit).

Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Disposal of used electric and electronic devices (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of these products.

Environment and health are endangered by a faulty disposal.

Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

Batteries

The included batteries are recyclable. Please dispose the batteries by using a dedicated battery container or return them to the specialist shop. Please dispose only discharged batteries to ensure environmental protection.

This device complies with the following EC directives:

- 2011/65/EC + 2015/863/EC RoHS directive
- 2012/19/EC WEEE directive
- 2014/53/EC RED directive

Compliance with the directives listed above is confirmed by the CE seal on the device.

CE compliance declarations are available on the Internet at www.humantechnik.com/service.

Technical specifications subject to change without prior notice.

Humantechnik's UK authorised representative:

Sarabec Ltd.,
15 High Force Road,
MIDDLESBROUGH TS2 1RH.
United Kingdom

Sarabec Ltd., hereby declares that this device complies with all UK statutory instruments.

UK declaration of conformity available from Sarabec Ltd.

Specifications

Power supply: 2 pcs 1.5 V AAA battery or
1.2 V rechargeable battery

Battery endurance:

approx. 800 h with 1.5 V alkaline battery

approx. 400 h with 1.5 V zinc carbon battery

approx. 300 h with 1.2 V rechargeable battery

Frequency: 868,425 MHz

Vibration: Eccentric motor, various
vibration sequences for identifying alerts

Height: 85 mm


Width: 60 mm

Depth: 28 mm

Weight: approx. 88 g

Operating range: 0° C – 40° C

Specifications subject to change without
prior notice.

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Humantechnik GHL AG