HUMANTECHNIK

signolux
Set

DE

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Seite 2

EN

User manual

Page 12

FR

Mode d'emploi

Page 22

NL

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Pagina 32

IT

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Pagina 42

ES

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Página 52

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We would like to congratulate you on purchasing a »signolux«. You have chosen a modern and reliable system. Kindly read the user manual carefully in order to start using the device properly and to familiarise yourself with all the possibilities of the system.

Standard delivery scope

Kindly check if all following parts are included:

- »signolux« receiver
- Doorbell- (Set 1) or call-pushbutton (Set 2)
- Battery (CR2032) for the doorbell pushbutton
- 4 batteries LR14
- 2 screw for the doorbell pushbutton
- Wall mount
- 2 screws for the »signolux« receiver
- 2 wall dowels
- Adhesive strips
- Operating manual
- Warranty card

Should any parts be missing, kindly contact your specialist dealer or the manufacturer directly.

Functioning principle

A »signolux« light signal system consists of at least one transmitter (e.g. doorbell pushbutton) and a »signolux« receiver. Up to 8 transmitters can be integrated.

If e.g. the doorbell pushbutton is pressed, a radio signal is transmitted to the »signolux« receiver.

The receiver then indicates the received signals acoustically and visually with sounds and light.

Safety instructions

- Prior to putting the device into operation and using it, kindly read this manual thoroughly and in full.
- Store this operating manual in a safe place so that it is also accessible for other users at any time.
- When using power tools always follow the instructions of the manufacturer and wear suitable protective equipment (e.g. protective goggles).
- Kindly check before starting to drill if there are any hidden electric cables and water pipes in the walls. In case of doubt we recommend using a cable detector.

Commissioning

1. Set 1: (A-2621-0) Inserting and replacing the battery of the doorbell pushbutton

In order to be able to insert the coin cell battery (type CR 2032) you first have to open the battery compartment.

- Place the doorbell pushbutton in such a way on the work surface that the two rotation axes of the transparent plastic face are at the rear (fig. 1 A and C).


Figure 1: Doorbell pushbutton

- b) Using a pointed object (e.g. biro) carefully press now the small pin on the left side (fig. 1 A).
- c) Now you should be able to push up and remove the transparent plastic face (fig. 1 A).
- d) You can now open the small screw (fig. 1 B) using a small Phillips screwdriver and remove the front of the doorbell pushbutton behind which there is the PCB with the battery.
- e) Insert the battery (coin cell CR 2032) into the plastic support on the PCB (fig. 2 D).


Figure 2: Battery compartment doorbell pushbutton

- f) If you would like to replace a battery, remove the old one with a small flat screwdriver from its plastic support by lifting the edge slightly (fig. 2 E). Afterwards you can insert a new battery. Kindly ensure correct polarity of the battery.
- g) Please refer to step 3 for the assembly of the doorbell pushbutton.

2. Label name tag and insert

- a) If you have not opened the housing yet, kindly proceed as described under step 1 a) to c).
- b) The transparent plastic face consists of two parts. The labelled name tag can be inserted between them (dimensions: 94 x 24 mm). Kindly ensure that the labelling is right way up. Kindly also keep in mind that the name tag has a cutout that does not cover the light diode or at least lets the light pass through.

- c) Put the two transparent parts of the plastic face back together so that the notch-side interlock.
- d) Please refer to step 3 for the assembly of the doorbell pushbutton.

3. Assembly of the doorbell pushbutton

- a) If you have changed only the name tag proceed to step e)
- b) Before assembling ensure that the small screw (fig. 1 B) can grip at the right spot. If required turn one part of the housing by 180°. Please make sure that the rubber gasket sits correctly within the housing.
- c) Insert the front face into the housing of the doorbell pushbutton and press both firmly together until the housing parts lock in place.
- d) Carefully drive in the small screw (fig. 1 B) with a small Phillips screwdriver.
- e) Now hold the doorbell pushbutton in such a way that the small screw (fig. 1 B) is on the right.
- f) At the right side of the pushbutton there is a small opening (fig. 1 C) where the hinge pin of the transparent plastic face can be inserted.
- g) You can now position the transparent plastic face in such a way until it has almost reached its final position. You only have to carefully press the small pin (fig. 1 A) and to lock the transparent plastic face in place.

4. Set 2: (A-2622-0) Inserting and replacing the battery of the call-pushbutton

In order to be able to insert the coin cell battery (type CR 2032) you first have to open the battery compartment.

- a) Remove the transparent front by lifting it off its clip on one side of the housing (fig. 3 A). If necessary, use a small flat-sided screwdriver to assist.


Figure 3: Call-pushbutton

- b) Now it should be possible to unhook the front from the the other side of the housing for removal.
- c) You can now open the small screw (fig. 3 B) using a small Phillips screwdriver and remove the front of the doorbell pushbutton behind which there is the PCB with the battery.
- d) Insert the battery (coin cell CR 2032) into the plastic support on the PCB (fig. 4 D).


Figure 4: Battery compartment call-pushbutton

- e) If you would like to replace a battery, remove the old one with a small flat screw-

driver from its plastic support by lifting the edge slightly (fig. 4 E). Afterwards you can insert a new battery. Please ensure correct polarity of the battery.

- f) Insert the front face into the housing of the call-pushbutton and press both firmly together until the housing parts lock in place. Please make sure that the rubber gasket sits correctly within the housing.
- g) Carefully drive in the small screw (fig. 1 B) with a small Phillips screwdriver.
- h) Hook the transparent plastic front onto the clip of one side, then press the other side down until the front completely locks onto the housing.

5. Inserting the batteries into the »signalux« receiver

- a) At the rear side there is a battery compartment the cover of which can carefully be pushed downwards so that the compartment opens (fig. 5).


Figure 5: Battery compartment of the
»signalux« receiver

- b) Now insert 4 alkaline batteries type LR14 (also called C batteries, included in the delivery scope). Ensure correct polarity of the batteries.
- c) Close again the battery cover by carefully inserting the cover and pushing it upwards.

Setting and programming


Figure 6: Receiver, upper side

6. Setting the signal mode

Using the signal mode slide switch (in fig. 6 "0 1 2") you can select the signal mode. Kindly keep in mind that the symbols are always displayed with coloured light diodes.

- switch to position 0:
signal tone, no flash,
- switch to position 1:
signal tones and flash
- switch to position 2:
no signal tones, flash only

7. Event icons

Figure 7: Event icons, control lamp and volume control

Each icon can be assigned to a transmitter (e.g. doorbell pushbutton) and can be redefined at any time.

The icons stand for:

Doorbell 1	Colour: Green
Doorbell 2	Colour: Green
Telephone	Colour: Orange

Person call	Colour: Yellow
Baby	Colour: White
Fire alarm	Colour: Red
Flood alarm	Colour: Blue
Alarm	Colour: Red

8. Setting the volume

The front side of the »signolux« receiver is protected by a screen that can be pushed upwards after loosening the lock on the top side of the »signolux«. Behind it you will find the rotating volume control with which you can set the volume of the signal tone.

Warning:

This device can generate sounds of up to 90 dB!

In order to protect your hearing, please set a low volume to begin with and only then adjust to the desired volume.

9. Programming signals and assignment to the icons

A transmitter (e.g. doorbell pushbutton) can be linked to any desired icon. This assignment is done by easy programming. Please keep in mind that a transmitter can only be programmed if it is not yet known in the memory of the »signolux«.

If required, you might first have to delete the transmitter from the memory (see step 11.)

- Keep the key »●« (fig. 6) pressed for some time. The symbols flash in sequence. Let go of the key as soon as you see the desired symbol.
- Trigger the signal now by using your transmitter (e.g. operate the doorbell pushbutton).
- Wait for about 20 seconds and test the function.

10. Setting the melody

You can adapt the melody for each icon:

- a) First trigger a signal (e.g. operate door-bell pushbutton).
- b) The icon assigned to the transmitter should illuminate and the previously chosen melody should be heard.
- c) Press and hold the melody setting button »♪« (fig. 6) for two seconds. Let go of the button and press it repeatedly until you have found the desired melody (you can choose from 16 melodies).
- d) Once the icon ceases to flash, the selected melody has been assigned.

11. Deleting transmitters from receiver memory

Do the following if you want to delete a transmitter from the memory:

- a) Keep the key »●« (fig. 6) pressed for some time. The icons light up in sequence. Let go of the key as soon as you see the desired icon.
- b) Now press key »●« and the melody setting button »♪« (fig. 6) at the same time. All icons light up briefly to confirm that the transmitter has been deleted from the receiver.

Please note that this will remove all assigned transmitters from the selected icon.

12. Battery level indication

The receiver informs about the state of the batteries of the receiver as well as the batteries of the transmitters:

- A sufficient battery level of the receiver is indicated by the red control lamp flashing every 5 to 6 seconds.
- A low battery level of the receiver is indicated by the red control lamp flashing in short sequence (once per second). Please replace the receiver batteries as soon as possible.
- A low battery level of a transmitter is indicated by the red control lamp shining together with the symbol of the according transmitter when it's signal is triggered. Please replace the batteries of the according transmitter as soon as possible.

Further functions

13. Checking the last signal, deleting the current signal

By pressing the »●« button (fig. 6) the last registered signal can be displayed again. You can also “delete” a signal that was just activated by pressing this button.

14. AUX-jack and power connector

The »AUX«-jack and the connector for the »signolux« power adapter (Item number A-2690-0, available as accessory) can be found on the right side of the receiver.

All Humanteknik modules can be connected to the AUX-jack (vibrating pillow etc, available as accessories)


Figure 8: Lateral connectors

a different mounting or installation location if the reception is not guaranteed.

15. Mounting the doorbell pushbutton

You can mount the doorbell pushbutton using the supplied strips of double-sided tape or the supplied set of screws.

16. Mounting the »signolux« receiver

The »signolux« receiver is suited either for free-standing installation or it can also be mounted using the supplied wall mount. To use it, fix the wall mount onto a wall (or rafter, cupboard etc.) When mounting onto a brick- or stone wall, you may need a power drill and stone bit to do so. Use the supplied wall plugs in this case. Once fixed, hang the »signolux« receiver onto the wall mount bracket.


Figure 9: Wall mount

Mounting

Important note: Kindly keep in mind that metal parts and other building materials may interfere with the transmission of the radio signal. As they can be either visible or hidden at the desired mounting location you should always carry out a functional test at the mounting location before starting the mounting process, and, if required, choose

Troubleshooting and solutions

Symptoms	Cause	Solutions
When programming a transmitter, a different icon lights up...	The transmitter you are trying to program is already known to the receiver.	See step 7. First the transmitter has to be deleted from the memory, before it can be assigned to a new symbol.
The red status LED light up with an alert...	This indicates a low battery in the according transmitter.	Insert a new battery in the transmitter.
The »signolux« receiver does not work...	Ensure that the batteries have been inserted correctly.	Ensure correct polarity of the batteries!
	Check if batteries type LR14 'C' have been inserted.	Use only alkaline batteries.
	It could be that the transmitter has not yet been programmed in the »signolux« receiver.	Programme it.
	It could be that the »signolux« receiver is out of reach of the transmitter.	Try using the transmitter or receiver at different location.
There is no ring tone...	switch position 2 (fig. 2)	Check the position of the switch (fig. 2) at the »signolux« receiver.
		Ensure that the volume is set to the lowest setting.
The flash light does not work...	switch position 0 (fig. 2)	Check the position of the switch (fig. 2) at the »signolux« receiver.

Symptoms	Cause	Solutions
The range is reduced...	Metal structures and other building materials can reduce the reach.	If possible do not mount the »signolux« transmitter or receiver near or at metal parts.
	Other devices may interfere with the radio signals, thus reducing the radio signals and the function of the »signolux« system.	
	Walls and ceilings reduce the range.	
	Low battery charge reduces the range.	Change the batteries every 12–18 months. In colder weather (below 5°C) you might have to change the batteries more frequently.

Maintenance and care

The »signolux« does not require any maintenance. If the unit does become dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents. Do not set up the unit where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Please make sure that the batteries are not exposed to excessive heat from such as sunlight, fire or anything similar.

Warranty

The »signolux« is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge. It is recommended that you send in the product in its original packaging, so keep the packaging for the duration of the warranty period.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

Batteries

The included batteries are recyclable. Please dispose the batteries through the battery container or the specialist shop. Please dispose only discharged batteries to ensure environmental protection..

Technical Data

	»signolux« receiver	doorbell- / callbutton
Work range	0 to 40°C	-10 to 40°C
Radio frequency	868,425 MHz	868,425 MHz
Range (open spaces)	200 m	200 m
Sound volume	90 dB at 1 m	-
HF capacity	-	<10 mW
Battery	LR14	CR2032
Battery service life (5 activations per day)	12 months	18 months
Power adapter connection	7.5 V DC / min. 800 mA	-
Dimensions (HxWxD)	144 x 114 x 41 mm	41 x 90 x 23 mm 80 x 40 x 16 mm
Weight	296 g	41 g

Hereby, Sarabec Ltd declares that this device complies with all relevant UK statutory Instruments.

The full declarations of conformity for this device are available from Sarabec Ltd., 15 High Force Road, MIDDLESBROUGH TS2 1RH. UK.

Sarabec's EU authorised representative:

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This device complies with the following EC standards:

- 2011 / 65 / EC + 2015 / 863 / EC RoHS-Directive
- 2012 / 19 / EC WEEE-Directive
- 2014 / 53 / EC RED-Directive

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at www.humantechnik.com.

Technical specifications subject to change without prior notice.

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