# Visual Signal System

# signolux

Doorbell-pushbutton

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We would like to congratulate you on purchasing a »signolux« signal system! You have chosen modern and reliable devices. Kindly read the user manual carefully in order to start using the device properly and to familiarise yourself with all the possibilities of the system.

### Package contents

Kindly check if all following parts are included:

- Doorbell-pushbutton
- Battery (CR2032)
- Double-sided adhesive strip
- User manual
- Warranty card

Should any parts be missing, kindly contact your specialist dealer or the manufacturer directly.

### Safety instructions

- Prior to putting the device into operation and using it, kindly read this manual thoroughly and in full.
- Store this operating manual in a safe place so that it is also accessible for other users at any time.
- When using power tools always follow the instructions of the manufacturer and wear suitable protective equipment (e.g. protective goggles).
- Kindly check before starting to drill if there are any hidden electric cables and water pipes in the walls. In case of doubt we recommend using a cable detector.

### Functioning principle

A »signolux« light signal system consists of at least one transmitter (e.g. doorbell pushbutton) and a »signolux« receiver. Up to 8 transmitters can be integrated. If e.g. the doorbell pushbutton is pressed, a radio signal (radio impulses on 418 MHz / 915 MHz) are transmitted to the »signolux« receiver. This one signals the received signals acoustically and visually with sounds and light signals.

## Setup

### 1. Inserting and replacing the battery

In order to be able to insert the coin cell battery (type CR 2032) you first have to open the battery compartment.

- Place the doorbell pushbutton in such a way on the work surface that the two rotation axes of the transparent plastic face are at the rear (fig. 1 A and C).


Figure 1: Doorbell-pushbutton

- b) Using a pointed object (e.g. pen) carefully press now the small pin on the left side (fig. 1 A).
- c) Now you should be able to push up and remove the transparent plastic face (fig. 1 A).
- d) You can now open the small screw (fig. 1 B) using a small Phillips screwdriver and remove the front of the doorbell pushbutton behind which there is the PCB with the battery.
- e) Insert the battery (coin cell CR 2032) into the plastic support on the PCB (fig. 2 D).


Figure 2: Battery compartment

- f) If you would like to replace a battery, remove the old one with a small flat screwdriver from its plastic support by lifting the edge slightly (fig. 2 E). Afterwards you can insert a new battery. Kindly ensure correct polarity of the battery.
- g) Please refer to step 3 for the assembly of the doorbell pushbutton.

## 2. Label name tag and insert

- a) If you have not opened the housing yet, kindly proceed as described under step 1 a) to c).
- b) The transparent plastic face consists of two parts. The labelled name tag can be inserted between them (dimensions: 94 x 24 mm). Kindly ensure that the labelling is right way up. Kindly also keep in mind that the name tag has a cutout that does not cover the light diode or at least lets the light pass through.
- c) Put the two transparent parts of the plastic face back together so that the notches interlock.
- d) Please refer to step 3 for the assembly of the doorbell pushbutton.

## 3. Assembly

- a) If you have changed only the name tag, proceed to step e)
- b) Before assembling ensure that the small screw (fig. 1 B) can grip at the right spot. If required turn one part of the housing by 180°. Please make sure that the rubber gasket sits correctly within the housing.
- c) Insert the front face into the housing of the doorbell pushbutton and press both firmly together until the housing parts lock in place.
- d) Carefully drive in the small screw (fig. 1 B) with a small Phillips screwdriver.
- e) Now hold the doorbell pushbutton in such a way that the small screw (fig. 1 B) is on the right.
- f) At the right side of the pushbutton there is a small opening (fig. 1 C) where the hinge pin of the transparent plastic face can be inserted.

- g) You can now position the transparent plastic face in such a way until it has almost reached its final position. You only have to carefully press the small pin (fig. 1 A) and to lock the transparent plastic face in place.

#### 4. Do not forget to register your doorbell-pushbutton with your »signolux«-receivers!

- a) Keep the function button of the according receiver pressed. The symbols light up in sequence. Release the button when the desired symbol is lit.
- b) Press the doorbell-pushbutton.
- c) Wait for approx. 20 seconds and test the function.

#### 5. Mounting the doorbell pushbutton

You can mount the doorbell pushbutton using the supplied strips of double-sided tape. Take care to adhere it to a smooth, dry and clean surface.

##### Important note:

Kindly keep in mind that metal parts and other building materials may interfere with the transmission of the radio signal. As they can be either visible or hidden at the desired mounting location you should always carry out a functional test at the mounting location before starting the mounting process, and, if required, choose a different mounting or installation location if the reception is not guaranteed.

#### Technical data

Work range	-10 bis 40°C
Radio frequency	418 MHz / 915 MHz
Range (open spaces)	up to 200 m
HF capacity	<10 mW
Battery	CR2032
Battery service life (5 activations per day)	18 months
Dimensions (HxWxD)	41 x 90 x 23 mm
Weight	41 g

This device complies to the following IC standards:  
47 CFR Part 15; 2013  
RSS-210 ISSUE 8

Technical specifications subject to change without prior notice.

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