Visual Signal System

signolux

Ruftaster
Oproepknop

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Deutsch 3

We would like to congratulate you on purchasing a »signalux« signal system! You have chosen modern and reliable devices. Kindly read the user manual carefully in order to start using the device properly and to familiarise yourself with all the possibilities of the system.

Package contents

Kindly check if all following parts are included:

- Person call-pushbutton
- Battery (CR2032)
- Double-sided adhesive strip
- User manual
- Warranty card

Should any parts be missing, kindly contact your specialist dealer or the manufacturer directly.

Safety instructions

- Prior to putting the device into operation and using it, kindly read this manual thoroughly and in full.
- Store this operating manual in a safe place so that it is also accessible for other users at any time.
- When using power tools always follow the instructions of the manufacturer and wear suitable protective equipment (e.g. protective goggles).
- Kindly check before starting to drill if there are any hidden electric cables and water pipes in the walls. In case of doubt we recommend using a a cable detector.

Functioning principle

A »signalux« light signal system consists of at least one transmitter (e.g. call-pushbutton) ad a »signalux« receiver. Up to 8 transmitters can be integrated.

If e.g. the doorbell pushbutton is pressed, a radio signal (radio impulses with 868,35 MHz) are transmitted to the »signalux« receiver. This one signals the received signals acoustically and visually with sounds and light signals.

Setup

1. Inserting and replacing the battery

In order to be able to insert the coin cell battery (type CR 2032) you first have to open the battery compartment.

- Remove the transparent front by lifting it off its clip on one side of the housing (fig. 1 A). If necessary, use a small flatsided screwdriver to assist.


Figure 1: Call-pushbutton

- b) Now it should be possible to unhook the front from the the other side of the housing for removal.
- c) You can now open the small screw (fig. 1 B) using a small Phillips screwdriver and remove the front of the call-push-button behind which the PCB with the battery. is located.
- d) Insert the battery (coin cell CR 2032) into the plastic support on the PCB (fig. 2 D).


Figure 2: Battery compartment

- e) If you would like to replace a battery, remove the old one with a small flat screwdriver from its plastic support by lifting the edge slightly (fig. 2 E). Afterwards you can insert a new battery. Kindly ensure correct polarity of the battery.
- f) Insert the front face into the housing of the call-push-button and press both firmly together until the housing parts lock in place. Please make sure that the rubber gasket sits correctly within the housing.
- g) Carefully set the small screw (fig. 1 B) with a small Phillips screwdriver.

- h) Hook the plastic front onto the clip of one side, then press the other side down until the front completely locks onto the housing.

2. Do not forget to register your call-pushbutton with your »signolux«-receivers!

- a) Keep the function button of the according receiver pressed. The symbols light up in sequence. Release the button when the desired symbol is lit.
- b) Press the call-pushbutton.
- c) Wait for approx. 20 seconds and test the function.

3. Mounting the call-pushbutton

You can mount the call-pushbutton, if needed, by using the supplied strip of double-sided tape. Take care to adhere it to a smooth, dry and clean surface.

Important note:

Kindly keep in mind that metal parts and other building materials may interfere with the transmission of the radio signal. As they can be either visible or hidden at the desired mounting location you should always carry out a functional test at the mounting location before starting the mounting process, and, if required, choose a different mounting or installation location if the reception is not guaranteed.

Technical data

Work range	-10 bis 40°C
Radio frequency	868,35 MHz
Range (open spaces)	up to 200 m
HF capacity	<10 mW
Battery	CR2032
Battery service life (5 activations per day)	18 months
Dimensions (HxWxD)	80 x 40 x 16 mm
Weight	41 g

This device complies to the following EC standards:

- 2011 / 65 / EC RoHS - Directive
- 2012 / 19 / EC WEEE - Directive
- 2014 / 53 / EC RED - Directive

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at

<http://www.humantechnik.com/service>

Technical specifications subject to change without prior notice.

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