HUMANTECHNIK

lisa
DS-1/RF

D

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Seite 2

Digitalwecker *lisa* DS-1/RF

GB

Operating Instructions

Page 10

Alarm clock *lisa* DS-1/RF

F

Mode d'emploi

Page 18

Réveil digital *lisa* DS-1/RF

NL

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Pagina 26

Digitale wekker *lisa* DS-1/RF

I

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Pagina 34



Sveglia digitale *lisa* DS-1/RF

E

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Página 42

Despertador *lisa* DS-1/RF

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

Congratulations on purchasing your alarm clock *lisa DS-1/RF*. We hope it will provide you with many years of service. You have chosen a modern and reliable system.

Please read the operating instructions through carefully to be able to start the unit correctly and to familiarise yourself with all of the systems features.

Standard components

Please check if all following components are included:

- *lisa DS-1/RF*
- 12 V power supply unit
- Operating instructions
- Warranty card

If any parts are missing please immediately contact your dealer or the manufacturer directly.

Functional principle

The alarm clock *lisa DS-1/RF* is an electronic clock with a digital time display. It wakes you up with an acoustic signal which is adjustable in audio frequency and volume.

The following accessories can be connected to the alarm clock:

- Vibrating pillow
- Flash module
- Switch module
- Acoustic module

Installation

Connect the power supply unit to the 12 V connector socket of the alarm clock and plug it into a 230 V power outlet. The display flashes as a sign that the time has not yet been set. The flashing disappears as soon as the time is set.

Time setting

Press  and **h** or **min** simultaneously to set the time.

Note

The time can be displayed in a 12 or 24 hour modus. The setting can be changed at the 12/24 h switch at the back of the alarm clock.

Setting the alarm time

The alarm clock has two alarm times that can be set separately.

Alarm time 1

Press  **1** and **h** or **min** simultaneously to set the desired alarm time (Duration of alarm signal 3 min).

Alarm time 2

Press  **2/Stop 24 h** and **h** or **min** simultaneously to set the desired alarm time. (Duration of alarm signal 3 min).

Note

After the installation only alarm time 1 is active.

Both alarm times can be activated separately. Press **1** and **2/Stop 24 h** simultaneously and push the snooze button. With every push of the snooze button one of the three possible settings is shown on the display.

You can choose between the following settings:

Alarm time 1 active:

Display indication: AL 1

Alarm time 2 active:


Display indication: AL 2

Alarm time 1 and 2 active:


Display indication: AL 1 2

In addition to the display indication the LED flashes as follows:

Alarm time 1 active:

Alarm time 2 active:

Alarm time 1 and alarm time 2 active:
LED lights continuously

Switching the alarm on/off

Activate/deactivate the alarm with the alarm on/off switch at the back of the *lisa DS-1/RF*.

If the switch is in the OFF-position, both alarm times are deactivated.

Setting of the audio frequency

The audio frequency of the alarm signal can be adjusted with the audio frequency control on the left side of the alarm clock. The audio frequency increases if the control is turned to the left and decreases if the control is turned to the right.

Note

The alarm tone is turned off if the audio frequency control is turned completely to the left.

Setting of the volume

The volume of the alarm signal can be changed with the volume control on the right side of the alarm clock. The volume decreases if the control is turned to the left, the volume increases if the control is turned to the right.

Snooze function

Pressing the snooze button interrupts the alarm. The alarm starts again after approx. 5 minutes.

Special features


- The alarm signal is repeated after 4 minutes (auto snooze) unless you press the snooze button or the Stop 24 h button.
- If a 9 V battery is inserted into the battery compartment the alarm clock operates even in case of a power failure. The display of the alarm clock turns off but the alarm signals are triggered.
- The inserted 9 V battery can be checked by pressing the snooze button. If "LO" is indicated on the display of the alarm clock, the battery has to be changed as otherwise the power failure back-up cannot be guaranteed.

The alarm clock *lisa DS-1/RF* is also a receiver for the *lisa*-RF-signalling system.

Setting the correct radio channel a special digital transmission process ensures that the radio signals are received error-free. This process also makes it possible to set 10 different radio channels. When supplied from the factory, the *lisa DS-1/RF* is set to channel 0. It is only necessary to set units to other channels if you have several different users operating units in the same building.

Set the channel you want by adjusting the channel selector switch on the back of the unit to the required channel number.

Use a small screwdriver to adjust this switch.

All transmitters and receivers within your system must be set to the same channel!

Transmission Signal Recognition:

The *lisa DS-1/RF* converts the various transmission signals into varying sequences of vibrating and signal displays.


Telephone

Telephone LED lights

Doorbell I

Doorbell LED lights

Doorbell II

Doorbell LED lights

Baby monitor

Babymonitor LED lights


Pager

Pager LED lights

Alarm

All LEDs light

Water alarm

Water alarm LED lights

Fire alarm

Fire alarm LED lights

P12 Pause 1,2 s

P14 Pause 1,4 s

P20 Pause 2 s

P30 Pause 3 s

P40 Pause 4 s

I Impulse

The indicator LEDs remain on for 40 seconds.

Touching the snooze button cancels any acoustic signal caused by *lisa* signals. The indicator LEDs light the normal time (40 seconds).

Recalling the most recent signal

The *lisa DS-1/RF* stores each signal it receives for a period of 3 minutes.

During this time you can call up the signal again by pressing the snooze button. The signal is then displayed on the indicator LEDs. An acoustic signal is not triggered.

Radio range

The transmitters are designed to have a radio signal range of 80 meters under optimum conditions. The following conditions may result in a reduced range:

- Signal transmission through building structures or vegetation
- Interference from TVs, computer monitors and mobile phones.
Under certain circumstances, this may prevent the receiver from functioning at all.
- Positioning of the transmitter at an unfavorable distance from surfaces which reflect radio waves, for example floors and walls.
This may lead to the signal wave being weakened or even cancelled out by the reflected wave.
- Metallic objects reduce the range due to the shielding effect of the metal.
- Particularly in towns and cities, there may be many other sources of radio waves which can disrupt the original signal.
- Having units operating at similar frequencies located close to one another may also result in mutual interference.

Maintenance and care

The *lisa DS-1/RF* does not require any maintenance. If the unit does become dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents.

Do not set up the unit where it will be exposed to full sunlight for long periods. In addition it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

Please make sure that the batteries are not exposed to excessive heat from such as sunlight, fire or anything similar.

Warranty

The *lisa DS-1/RF* is a very reliable product. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer.

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by

people who are not authorised to do so (destruction of the seal on the unit).

Repairs under warranty are only carried out providing you have filled in and returned the enclosed warranty card from the dealer.

Always specify the product number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system). The symbol on the product or the packa-

ging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

Technical data

Power supply: 12 V power supply unit at 230 V, 50 Hz or 9 V battery
Power consumption: approx. 2.5 Watt


This device satisfies the following EU directives:

- 2011/65/EG RoHS - directive
- 2002/96/EG WEEE - directive
- 1999/5/EG R&TTE - directive
- EN 301 489-3 V1.6.1
- EN 60950-1: 2006 + A11: 2006 + A1: 2010 + A12: 2011
- EN 55022: 2010
- EN 300 220-2 V2.4.1

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at www.humantechnik.com.

Technical specifications subject to change without notice.

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