DE

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Seite 2

Ringschleifenverstärker

EN

Operation Instructions

Page 10

Loop amplifier

FR

Mode d'emploi

Page 18

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
1. Overview of the connections and operating elements

1.1 Front

- 1 Input signal indicator
- 2 Loop indicator
- 3 Regulator to set the corresponding input sensitivity
- 4 Treble
- 5 Regulator to set the loop power
- 6 Power supply indicator
- 7 Earphone jack
- 8 Volume control for the earphone jack

1.2 Rear

- 9 Power supply 230 V / 50 Hz
- 10 Screw terminal for the connection of the loop
- 11 Line OUT – Monitor jack
- 12 Line IN - Input for devices with Line level (25 mV - 4 V)
- 13 Line IN - Inputs with symmetrical XLR jacks. Per jumper: variable input sensitivity between 0.5 mV - 100 mV for microphones and 25 mV - 4 V for devices with Line level. Phantom voltage can be activated.

We would like to congratulate you on the purchase of your *PROLOOP C* loop amplifier. You have chosen a product that skillfully combines appealing design with high technical performance.

This handbook describes exactly how the amplifier works and how it is to be used.

Please read through the operating instructions carefully before using the amplifier.

We also reserve the right to make changes to incorporate new product developments at any time.

1.3 Package contents

The following items are included:

- *PROLOOP C* induction loop amplifier
- Power cord
- 2 Mounting plates
- Operating instructions
- Guarantee card

Should any of these items be missing, please contact your local retailer or the manufacturer.

2. Notes on safe operation

This device complies with all the necessary EU directives and bears the CE seal of approval. To make sure that your device enjoys a long service life, please pay attention to the following information.

Important notes

- **To prevent accidents or personal injury caused by electrical shocks, never place any type of container filled with a liquid, such as a vase, on the device.**
- **Never pull on the power cord to remove the plug from the wall outlet; always pull the plug.**
- **Do not operate the device near heat sources or in rooms with high humidity. (Operating temperature range 0-40°C).**
- **Do not cover the air vents so that any heat generated by the device can be dissipated by air circulation.**
- **Depending on how it is used, the device can develop high temperatures in the heat sink and at the air vents. Be careful not to touch these elements – risk of burns !!!**

Important notes

Do not use the device:

- if the power cord or the plug are defective or damaged.
 - if the device is not functioning properly or is defective.
 - if the device has fallen or been damaged in any other way.
-
- Always make sure to pull the plug before attempting to clean the device.
 - To clean the device, simply use a soft, dry cloth – never use any kinds of chemicals or water.
 - Should a repair become necessary, have this done by an authorized technician only.

3. Applications

The *PROLOOP C* amplifier was designed for use in professional venues where high reliability and safe operation are essential. The 100 % short-circuit-proof amplifier with programmable, symmetrical XLR inputs and extremely stable output power meets these requirements. The AGC function attenuates the signal level of such constant signals as oscillation and sine curves to - 10 dB, thus providing for consistent field strength and stable sound.

Inductive listening systems are used for different applications, e.g. for assisted listening for the hearing impaired in churches, theaters, cinemas, waiting rooms and lounges, as interpreting systems, for lectures in museums, exhibitions, etc..

3.1 How it works

The induction loop amplifier generates an alternating electromagnetic field around the loop cable. The loop consists of an insulated wire, which is usually laid singly (or doubly) along the side of the room.

When the hearing-aid wearer set the hearing aid to the “T position”, voltage is induced into the T coil. This voltage depends on the amplifier’s input signal. Now the signals are amplified by the hearing aid and converted into sound. The alternating electromagnetic field inside the loop allows the hearing-aid wearer to move around freely in the room and to receive the voice or music signals as high-quality, pleasant sound.

3.2 PROLOOP C loop amplifier performance features

- For areas up to 170 m²
- 4 A RMS output power
- Short-circuit-proof
- Automatic reset for the fuse
- 2 programmable, symmetrical XLR inputs and 1 RCA Line IN input
- Line IN or microphone IN with phantom voltage can be adjusted separately for Line IN 1 and Line IN 2
- High safety thanks to the output voltage and power controlled AGC, which reduces problems with acoustic feedback
- Headphones connector (front side) and Line OUT RCA connector (rear side)
- Treble control to compensate for treble loss due to reinforcements in walls and floors

4. Installation

This device operates with a dangerous voltage of 230 V. Never attempt to open the device yourself and never insert anything into the air vents! You could receive a severe electrical shock!


Advice

1. Before beginning to install the device, you should be sure to check the area for any magnetic interference. Such interference can have a negative affect on the operation of the device or make it impossible altogether. Interference can be caused by transformers, power lines, reinforced concrete or floor heating systems.
2. If the loop is to be installed in a pipe, then make sure that the pipe is made of plastic.
3. Do not lay any input signal cables next-to or along the loop wire.
4. Avoid using dynamic microphones inside the loop. Always chose symmetrical electret or capacitor microphones to prevent acoustic feedback.
5. If it is not possible to lay the loop in a circle or in a "figure 8" (double circle), then a special loop design calculated by an experienced technician will be necessary. (Please do not lay any comb or meander structures!)

4.1 Connecting the induction loop


Connect the ends of the loop wires to the terminals (10) on the back of the amplifier.


Please note the following information:

- The loop amplifier must be located outside of the loop.
- Twist the cable section between the amplifier and the loop (feed cable to the loop) to prevent interference.
- Use a suitable measurement device (e.g. ohm meter) to make sure that the loop doesn't have a short-to-ground (no damage on the insulation of the loop wiring)!
- The distance from the loop level to the normal listening level should be approx. 1.20 m to 2.00 m.

4.2 Notice for 19" rack installation

Remove the 3 screws on the forward area of the left and right sides of the *PROLOOP C*.

Use the removed screws to fix the 19" rack mounting plates on the forward areas of the left and right side of the *PROLOOP C*.

5.1 Configuring the amplifier

Two audio sources can be connected to the symmetrical XLR jacks (13). If a microphone is connected, the corresponding jumper must be configured accordingly.

If phantom-fed microphones are used, a 40 V phantom voltage can be added.

Input IN3 (12) is exclusively designed for devices with Line level and cannot be configured with the jumpers.

In the standard configuration, input IN1 is for “microphone” with phantom voltage, input IN2 for “Line IN” without phantom voltage.

If you want to alter the standard setting, **first pull the plug out of the wall outlet** and then remove the device cover and configure the device according to the following directions:


Figure: Jumper positions

Jumper	Function	Open	Closed	Default setting
IN1 S1 / S2 S3	Phantom voltage IN1 * Sensitivity IN1	Off Microphone	On Line In	On (closed) Microphone (open)
IN2 S4 S5 / S6	Sensitivity IN2 Phantom voltage IN2 *	Microphone Off	Line In On	Line In (closed) Off (open)
IN3 -	-	-	-	-

* With phantom voltage switched on, **do not** connect a microphone or a Line-IN device with an unbalanced jack - this could otherwise result in damage to the device!

6.1 First-time operation

Before you switch on the device, please turn all the control knobs as far to the left as possible.

Please note that the amplifier does not have an on/off switch and will switch on as soon as the power cord is plugged into the wall outlet.

Now connect the power cord that came with the device. The indicator for the power supply (6) will light up.

Before you adjust the Line level, please don't forget to send a test signal to the inputs.

Adjust the level of the Line-In inputs one after the other so that the input signal control lamp (1) lights up on the front for the program peaks of the signal source.

Then increase the power control (5) until the loop control lamp (2) lights up on the front of the device.

Use the earphone jack (7) to listen directly and, if necessary, adjust the treble with the treble control (4).

6.2 Adjusting the field strength

You will need a field strength meter (special accessory, e.g. „PROLOOP FSMplus“, Item No. A-4292-0) to calibrate the system correctly.

7. Warranty

The induction loop amplifier exhibits very high operating safety. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly. This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

8. Disposal

of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units. You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal. Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

9. PROLOOP C specifications

Input

Connectors	2 x XLR-jacks balanced 1 x RCA
max. input level	setting by jumper, variable sensitivity level (0,5 mV - 100 mV) Mikrofon / (25 mV - 4 V) Line
Phantom voltage	activation by jumper setting

Line Output

Connector	1 x RCA
-----------	---------

Power amplifier output

Connection	Screw terminals for induction loop
max. loop power	4 A RMS
max. loop voltage	31 V
Coverage	170 m ² acc. to IEC 60118-4 single loop, free field

Transmission gain

Frequency range	100-5000 Hz (± 3 dB)
Distortion	Distortion < 1%

Power supply

Performance	230 V AC / 50-60 Hz
-------------	---------------------

Indicators

Power supply	1 green LED
Power induction loop	1 green LED
Input signal	1 green LED

Dimensions W x H x D

295 x 64 x 205 mm

Weight

approx. 3,6 kg

Colour

black

This device complies to the following EC directives:

- 2011/65/EC RoHS-directive
- 2002/96/EC WEEE-directive
- 2006/95/EC Low voltage directive
EN 60065: 2002 + A1: 2006 + A11: 2008 + A2: 2010 + A12: 2011
- 2004/108/EC Electromagnetic Compatibility
EN 55013: 2013 • EN 55020: 2007 + A11: 2011
EN 61000-3-2: 2006 + A1: 2009 + A2: 2009 • EN 61000-3-3: 2013

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at www.humantechnik.com.

Technical specifications subject to change without notice.

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Sarabec Ltd