**DE**

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Seite 2

Ringschleifenverstärker

**EN**

**User manual**

Page 10

Loop amplifier

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
We would like to congratulate you on purchasing your new *PROLOOP LOS* induction loop amplifier. You have chosen a product that skillfully combines appealing design with high technical performance.

This handbook describes exactly how the amplifier works and how it is to be used. **Please read this user manual thoroughly** before using the amplifier. We also reserve the right to make changes to incorporate new product developments at any time.

## 1. Connections and operating elements

### 1.1 Front panel

- (1) Power supply indicator
- (2) Input sensitivity indicator
- (3) Regulator to set required input sensitivity
- (4) Regulator to set slave loop current
- (5) Volume control for headphone jack
- (6) Slave headphone jack
- (7) Level indicator for the loop current (slave)
- (8) Regulator to adjust the master loop current
- (9) Volume control for headphone jack
- (10) Master headphone output jack
- (11) Level indicator for loop current (master)

### 1.2 Rear panel

- (12) Power supply (230V / 50Hz)
- (13) Screw terminals for the connection of the loop (master/slave)
- (14) LINE OUT, monitor output (master/slave)
- (15) 4 pole DIP switch
- (16) IN3 RCA LINE IN input jack
- (17) IN2 RCA LINE IN input jack
- (18) IN1 LINE IN jack Combo XLR (symmetrical) / 6.3mm cinch (asymmetrical)

### 1.3 Package contents

Please make sure that all the following parts are in the box:

- *PROLOOP LOS* induction loop amplifier
- Power cord
- Operating instructions
- “Induction Listening” sticker
- Guarantee card

If anything is missing, please contact your retailer or the manufacturer.

### 1.4. General information

The *PROLOOP LOS* amplifier has been designed for use in professional venues where high reliability and safe operation are essential. The 100% short-circuit-proof amplifier with its extremely stable output power meets these requirements.

### 1.5. Notes on safe operation

This device complies with all the necessary EU directives and bears the CE seal of approval. To make sure that your device enjoys a long service life, please observe the following information.

#### **Important Information**

To prevent accidents or personal injury caused by electrical shocks, never place any type of container filled with a liquid, such as a vase, on the device.

Never pull on the power cord to remove the plug from the wall outlet; always grasp the plug directly and pull it out.

Do not operate the device near heat sources or in rooms with high humidity (operating temperature range 0-40°C).

Do not cover the air vents so that any heat generated by the device can be dissipated by air circulation.

Depending on how it is used, the device can develop high temperatures in the heat sink and at the air vents. Be careful not to touch these elements by accident – risk of burns !!!

Do not use the device:

- if the power cord or the plug are defective or damaged.
- if the device is not functioning properly or is defective.
- if the device has been dropped or damaged in any other way.

To clean the device, simply wipe it with a soft, dry cloth – never use any kinds of chemicals or water.

Should a repair become necessary, have this done by an authorized technician only.

## 2. Why LOS?

Standard induction loop systems create a magnetic field that usually goes beyond the »key supply area«. This horizontal and vertical »overspill« can still be audible in areas up to three or four times outside the actual induction loop, which tends to limit the applications for which this technology can be used. Low-overspill systems (LOS), designed and configured by qualified experts, significantly reduce overspill and make it possible to use induction loop systems in directly adjacent rooms.

The PROLOOP LOS is a complete system for the operation of low-overspill systems. It includes all the electronics required in a single housing.

### 2.1 LOS induction loops

The ideal design and configurations are always dependent on the venue details, i.e. the floor plans and the sizes of the rooms serviced, the number of internal induction loops and the demands made on low overspill.

To determine the corresponding values for the »tailored« low-overspill layout and the optimal level adjustments in each case, it is necessary to use complicated calculation and installation methods. That's why we recommend that you use the services provided by HUMANTECHNIK. We support architects, technicians and installation experts qualified in the planning of LOS installations.

## 3. Installation notes

This device operates with a dangerous voltage of 230 V. Never attempt to open the device yourself and never insert anything into the air vents! You could receive a severe electrical shock!

### IMPORTANT INFORMATION

When installing the induction loop amplifier, make sure that the air vents on the front and rear panels are not covered or blocked in any way.

It is normal for the amplifier to become warm while it is in operation. It is important that the ventilation of the device is not impeded in any way. If the incoming air supply is limited for some reason, it will be necessary to supply forced cooling.

Make sure that the amplifier is securely mounted and cannot slip or fall.

To prevent damage to the inner components, do not insert any screws or drill any holes into the device housing.

Make sure that the voltage indicated on the rear panel of the device is the same as your local voltage.

All connections should be made by a qualified technician only and only when the amplifier is switched off.

### 3.1 Installation of the induction loop

To use the LOS amplifier, the induction loop must be installed in a special way. In detail, two separate loops are used (master/slave), which consist of eight-piece segments lined up one after the other. The exact proportions and the exact number of segments depend on the demands made on low overspill as well as the size of the room. The smaller the segment, the smaller the overspill will be.

The master loop is laid over the entire area to be serviced. The slave loop is always one segment shorter than the master loop and doesn't begin until after the first segment half of the master loop. The individual segment lengths should correspond to the distances between the seating rows.


We recommend using a single-strand cable, e.g. H07V-K with a cross section of 4mm<sup>2</sup>, or the copper ribbon cable (A-4937-0).

Advantages of a LOS induction loop design


- No loss of field strength in the centre of the induction loop
- No horizontal overspill
- No signal holes inside the loop

### 3.2 Example of a LOS loop measuring 16x19 metres, approx. 300m<sup>2</sup>

Master-Loop


Slave - Loop

*Fig.1 Example of a LOS loop*

#### 4. Connecting the induction loop

- The loop amplifier has to be located outside of the loop.
- Twist the section of cable between the amplifier and the loop (supply to the induction loop) to reduce interference. (See Fig.2)
- Use a suitable measuring instrument (e.g. ohmmeter) to make sure that the loop is not grounded (no defect in the insulation of the loop system). Typical resistance values for a 300m<sup>2</sup> loop lie between 0.5 and 2 ohms. Since the master loop is longer, it also has a greater resistance.
- Connect the cable ends of the loop to the corresponding outer terminals (13) on the rear panel of the amplifier. Make sure that you do not mix up the ends of the master and the slave cable.

*Fig.2 Connecting the induction loop*

#### 5. Configuring the audio inputs

Up to three audio sources with line level can be connected to the audio inputs (15,16,17). The input IN3 does not have a sensitivity regulator on the front panel because it has a predefined level. Any devices equipped with a line level can be connected, e.g. mixing desks, tuners, CD and MP3 players, TVs, receivers, DVD players or video recorders.

If you want to connect a microphone, please see 5.1 Connecting a microphone – Step by step

## 5.1 Connecting a microphone – Step by step

To connect a microphone, you have to use the IN1 input (see Fig. 3). The Combo jack allows you to plug in a microphone with an XLR plug (symmetrical) or with a 6.3mm cinch plug (asymmetrical).



Abb.3 Connect a microphone to IN1 input


1. First disconnect the plug from the wall outlet.
2. Slide the DIP switch 1 and 2 to position „OFF“ (see 5.2 Overview, Microphone level)
3. If you use a microphone with phantom voltage, also slide the DIP switch 3 to position „ON“ (see 5.2 Overview, Microphone level with phantom voltage).  
**Important:** When the phantom voltage is switched on, do not connect a microphone with an asymmetrical output or a LINE IN device, as this can result in damage to these devices.
4. Plug the power cord back into the wall outlet.

## 5.2 Overview – Setting options for the IN1 jack

Line level:


DIP switch 1 and 2 at position „ON“ - Line level  
**Important:** DIP switch 3 at position „OFF“ - phantom voltage switched off  
DIP switch 4 at position „ON“ - AGC switched on


DIP switch 1 and 2 at position „ON“ - Line-Pegel  
**Important:** DIP switch 3 at position „OFF“ - phantom voltage switched off  
DIP switch 4 at position „OFF“ - AGC switched off

Microphone level with phantom voltage:


DIP switch 1 and 2 at position „OFF“ - Microphone level  
DIP switch 3 at position „ON“ - phantom voltage switched on  
DIP switch 4 at position „ON“ - AGC switched on


DIP switch 1 and 2 at position „OFF“ - Microphone level  
DIP switch 3 at position „ON“ - phantom voltage switched on  
DIP switch 4 at position „OFF“ - AGC switched off

Microphone level without phantom voltage:

DIP switch 1 and 2 auf position „OFF“ - Microphone level  
DIP switch 3 at position „OFF“ - phantom voltage switched off  
DIP switch 4 at position „ON“ - AGC switched on

DIP switch 1 and 2 at position „OFF“ - Microphone level  
DIP switch 3 at position „OFF“ - phantom voltage switched off  
DIP switch 4 at position „OFF“ - AGC switched off

## 5.3 Switching on and off the AGC function (Automatic Gain Control)

To switch on the AGC function, slide the DIP switch 4 up to position „ON“. The AGC function will be switch off, if you slide the DIP switch 4 back down.

## 6. Adjusting the amplifier

In order to adjust a LOS system, you will require a field strength gauge, e.g. the FSM (A-4291-0).

1. Before you switch on the amplifier, turn all the dials as far to the left as possible. Please note that the amplifier does not have an ON/OFF button and is switched ON as soon as it is plugged into the wall outlet.
2. Plug in the power cord supplied with the amplifier. The power supply indicator (1) should now light up.
3. Before you adjust the line level, don't forget to apply a test signal to the corresponding input. We recommend using our test CD, which comes with the FSM (A-4291).  
Use a 1KHz sine signal (track 3 on the CD) and set the volume control (master fader) on the mixing desk or the audio amplifier as it is normally used at your events.
4. Adjust the input sensitivity regulator (3) of the LINE IN input used so that the control indicator (2) on the front panel of the device lights up constantly.
5. Turn the master power regulator (8) to the central setting. The control indicator (11) on the front panel should now light up.
6. Set the FSM to "Norm" and "0dB scale" and hold it vertically to the induction loop level. The ideal height depends on what the induction loop system is being used for; typically, this might be 1.20m (for seated audiences) or 1.70m (for standing audiences).
7. Check the average field strength by running diagonally over the induction loop area and determine the average of the values indicated. The average value has to measure -12dB ( $\pm 3$ dB).  
If you do not obtain this value, adjust the corresponding control until this value is reached. Please note that every time you cross a cable, the field strength will drop briefly. This is normal as long as only one loop is in operation.
8. Once the master loop is adjusted, unplug the amplifier from the power supply and disconnect the master loop. Important: Do not alter any of the controls mentioned in the above.
9. Plug the amplifier back into the power supply and repeat Steps 5 to 7 for the slave loop or the slave loop current regulator (4).
10. Once the slave loop is adjusted, reconnect the master loop.
11. For the final check, use a normal audio signal, such as music or voice via the microphone and check the field strength with the FSM. The value should read 0dB ( $\pm 3$ dB) in the signal peaks.  
If you do not read these values, repeat the measurement process starting from Step 3 and set both controls higher than for the first test.
12. If the values lie within the stipulated range, the amplifier is set to the ideal levels. Use the headphone jack of the FSM to evaluate the sound of the signal with a pair of headphones.

## 7. Warranty

The *PROLOOP LOS induction loop amplifier* exhibits very high operating safety. If any malfunctions should occur, despite correct operation, please contact your retailer or the manufacturer. The guarantee covers repairs and shipping free of charge. The device must be returned in its original packaging – do not discard it. The guarantee becomes invalid if the device is operated incorrectly or any repairs are attempted by unauthorized persons (device seal is broken). Repairs will only be made on guarantee if you also include the correctly completed guarantee card and a copy of the retailer's invoice/receipt. **The device number must always be included.**

## 8. Environmental protection / disposal

Disposal of used electrical and electronic devices (applicable in EU countries and other European countries with a separate collection system for these devices.) The symbol on the product or its packaging indicates that this product should not be thrown into the regular household waste; it should be taken to a recycling centre for electrical and electronic devices. When you dispose of this product correctly, you make a contribution towards protecting the environment and your health. Incorrect disposal can pose a threat to the environment and your health. Material recycling also helps to reduce the use of raw materials. Further information on recycling this product is available from your community, the municipal waste disposal department or the retail outlet at which you bought the product.

## 9. PROLOOP LOS specifications

|                                |   |
|--------------------------------|---|
| Power supply                   | 230 – 240 V AC, 50Hz, 25 – 700 W  |
| Coverage                       | 300 m <sup>2</sup> acc. to IEC 60118-4  |
| <b>Induction loop output</b>   |   |
| Max. power                     | 60 A from peak to peak, 1-5ms, 1 kHz, 2 x 13 A RMS, continuous 1 kHz, short-circuit-proof   |
| Max. voltage                   | 47 V from peak to peak  |
| Frequency range                | 100 – 5000 Hz ( $\pm 3$ dB)   |
| Distortion                     | < 1 %   |
| Cable connection               | 2 screw connections on the rear panel of the device   |
| <b>Outputs</b>                 |   |
| Master LINE OUT                | 0 dBm, RCA (with AGC function)  |
| Slave LINE OUT                 | 0 dBm, RCA (with AGC function)  |
| <b>Inputs</b>                  |   |
| IN 1                           | Kombinierter XLR / 6,3 mm Verbindungsstecker<br>LINE-Empfindlichkeit 37 mV – 10V RMS / 10kOhm<br>MIC-Empfindlichkeit 2,5 mV – 10V RMS / 10kOhm,<br>mit zuschaltbarer Phantomspannung 9 – 20 V |
| IN 2                           | Dual RCA connectors, 45 mV – 10V RMS / 10kOhm   |
| IN 3                           | Dual RCA connectors, 45 mV – 10V RMS / 10kOhm   |
| <b>Controls and regulators</b> |   |
| Induction loop regulators      | 0 – 300m <sup>2</sup> (one potentiometer for the master and slave output each)  |
| Input Level                    | IN 1 and IN 2 can be adjusted separately  |
| <b>Indicators</b>              |   |
| Power supply                   | 1 green LED   |
| Input level                    | 1 green LED   |
| Induction loop current         | 2 green LEDs  |
| Induction loop monitor         | Two 6.3-mm jacks to connect headphones  |
| Dimensions                     | 88 x 438 x 280 mm (H x W x D)   |
| Colour                         | Black   |
| Weight                         | 9,3 kg  |

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For other service-partners  
in Europe please contact:

**Humantechnik Germany**

Tel.: +49 (0) 76 21/ 9 56 89-0  
Fax: +49 (0) 76 21/ 9 56 89-70  
Internet: [www.humantechnik.com](http://www.humantechnik.com)  
E-mail: [info@humantechnik.com](mailto:info@humantechnik.com)

# AUDIOropa

... a Division of the Humantechnik Group

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**Sarabec Ltd**