DE

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Class D Ringschleifenverstärker

EN

User manual

Class D loop driver

INHALT

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INDEX

1. Introduction.....	17
2. PRO LOOP NX7 LOS.....	17
2.1 Description.....	17
2.2 Performance range.....	17
2.3 Contents of package.....	17
2.4 Advice and safety.....	18
3. Function.....	18
4. Indicators, connectors and controls.....	18
4.1 Indicators.....	18
4.2 Loop output.....	18
4.3 Front panel and controls.....	19
4.4 Rear panel and connectors.....	20
4.5 Rack installation.....	21
4.6 Adjustments and connectors.....	22
Maintenance and care.....	24
Warranty.....	24
Disposal.....	24
Specifications.....	25

1. Introduction

Thank you for having purchased the »PRO LOOP NX7 LOS« Class D loop driver! Please take a few moments to read this manual. It will ensure you the best use of the product and many years of service.

2. PRO LOOP NX7 LOS

2.1 Description

The PRO LOOP NX series consists of Class D loop drivers made to equip rooms with audio support for people with hearing loss.

2.2 Performance range

The »PRO LOOP NX7 LOS« belongs to a generation of induction loop drivers with high performance and efficiency. With this device it is possible to establish installations according to the international standard IEC 60118-4.

2.3 Contents of package

Please check if the following pieces are included in the package:

- PRO LOOP NX7 LOS induction loop driver
- Power cable 1.5 m, connectors CEE 7/7 - C13
- 2 pieces 3-point Euroblock-connectors
- 1 piece 2-point Euroblock-connectors, input
- 2 pieces Neutrik NL4FX speakON cable socket for loop output
- Mounting screws
- Adhesive loop-indication signs

Should any of these items be missing, please contact your retailer.

2.4 Advice and safety

- Never pull on the power cord to remove the plug from the wall outlet; always pull the plug.
- Do not operate the device near heat sources or in rooms with high humidity.
- Do not cover the air vents so that any heat generated by the device can be dissipated by air circulation.
- An installation must be carried out by qualified personnel.
- The device must be out of reach of unauthorized persons.
- The device is only to be used for operating inductive loop systems.
- Install the device and its wiring in such a way that there is no danger, e.g. by falling or tripping.
- Connect the loop driver only to wiring which complies to IEC 60364.

3. Function

An inductive listening system basically consists of a copper wire connected to a loop amplifier. Connected to an audio source, the loop amplifier generates a magnetic field in the copper conductor. The listener's hearing aids receive these inductive audio signals wirelessly in real time and directly in the ear - free from distracting ambient noise.

4. Indicators, connectors and controls


4.1 Indicators

The function status of the loop amplifier is continuously monitored. If an error occurs within the amplifier unit, the loop conductor or the inputs, this is indicated by corresponding LEDs on the front panel. The error status is also given out via the fault relay at the same time.

4.2 Loop output


The »PRO LOOP NX7 LOS« is constructed with two channels and thus has two loop outputs.

4.3 Front panel and controls

- 1 Power-LED - Indicates operation
- 2 IN 1: For adjusting the Mic/Line level of input 1
- 3 IN 2: For adjusting the Mic/Line level of input 2
- 4 IN 3: For adjusting the 100V level of input 3
- 5 Compression: Display of the level reduction in dB, in relation to the input signal
- 6 MLC (Metal Loss Correction) Compensation of frequency response due to metal influence in the building
- 7 Overload indicator »Protect / Clip«-LED (red) for Loop A
»Loop«-LED (yellow) Error indicator for Loop A
- 8 Loop current regulator Loop A
- 9 Output level display for Loop A and Loop B
- 10 Overload indicator »Protect / Clip«-LED (red) for Loop B
»Loop«-LED (yellow) Error indicator for Loop B
- 11 Loop current regulator Loop B
- 12 Loop A or Loop B selector switch for headphone output
- 13 Headphone output for monitoring the outgoing audio signal

4.4 Rear panel and connectors

14 INPUT 1 Combo: Mic or Line input

15 INPUT 1 selector switch: Switches between microphone- or line level.

16 Phantom power: Activates 48V phantom power on INPUT 1

17 INPUT 2: Mic or line input

18 Input selector switch for INPUT 2. Switches between „Mic“ or „Line“ input types.

19 Phantom power switch. Activates the supply current for microphones on INPUT 2.

20 Link input: Used in combination with the output of another PRO LOOP NX loop driver

21 Link output: Used in combination with the input of another PRO LOOP NX loop driver

22 Link output 90°: Used in combination with the link input of another PRO LOOP NX driver

23 0°/90° switch: Switches a 90° phase shift on or off in the Loop B output

24 100V Input: Input for audio signals of a 100V system

25 Fault: Output relay for the fault status of the device as a potential free circuit.

26 Loop A: NL4 Speaker Twist connection of loop output channel A


27 Loop B: NL4 Speaker Twist connection of loop output channel B

28 Mains socket with fuse holder

29 Mains switch for switching the unit on and off.

4.5 Rack installation

Leave a space of 1U (44mm) above the loop driver to ensure air circulation for its cooling.

Attention, Warning, Danger :

The loop driver contains an advanced protection circuit, which allows it to reduce the power output to maintain safe operating temperatures.

Insufficient ventilation may cause the amplifier output power to be reduced during normal operation (indicated by the red PROTECT LEDs lighting up).

To reduce the risk of thermal limitation and to allow proper heat dissipation, it is recommended to keep the space directly above and behind the device clear.

4.6 Adjustments and connectors

4.6.1 Loop connectors

The induction loops are connected via the NL4 Speaker Twist connectors for loop A and loop B located on the rear panel of the driver.

These have four poles, with 1+ parallel to 2+ and 1- parallel to 2-.

4.6.2 Audio inputs

Audio sources connect via the 3 inputs of the driver provided for this purpose.

The driver has 3 types of input :

INPUT 1: Line or microphone level

INPUT 2: Line or microphone level

INPUT 3: 100 V level

4.6.3 100V priority input

INPUT 3 (100V) on the rear panel of the driver is designed to prioritize 100 V PA systems. If several active audio sources are fed to the inputs of the driver, it will always be INPUT 3 that has priority over the others.

Incoming signals on INPUT 3 will mute INPUT 1 and INPUT 2.

4.6.4 Link connectors

On the rear of the unit you will find 2 sockets for 6.35mm jacks. With these balanced connections, several loop drivers can be connected with each other.

4.6.5 Line output

It is possible to connect a recording device and other gear to the line output of the loop driver. The line output of the amplifier is the 6.35mm jack marked »OUT«, which is also used to link several loop drivers when installing complex systems.

A further 3.5 mm jack socket can be found on the front of the device, for easy monitoring of the audio output via headphones.

4.6.6 Power supply

PRO LOOP NX drivers use a direct power supply of 100 - 265 V AC - 50/60 Hz.

4.6.7 Fault status relay


The error status of the loop driver can be read out via the potential-free relay output.

4.6.8 Terminal assignment:

All connectors are electronically balanced..


XLR

- Pin 1 : Ground/Shield
- Pin 2 : not inverted / hot (+)
- Pin 3 : inverted / cold (-)

6,35 mm Stereo jack

- Tip not inverted / hot (+)
- Ring inverted / cold (-)
- Shaft Ground/Shield

4.6.9 Power on / off

The unit is powered up using the switch on the rear panel of the loop driver. If the device is powered, the Power LED on the front panel lights up.

To switch the unit off, press the power switch again.

Please note that the device enters a stand-by mode when switched off.

Disconnect the mains plug from the socket in order to switch off the device completely.

Maintenance and care

The »PRO LOOP NX7 LOS« does not require any maintenance under normal circumstances. If the unit becomes dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents. Do not place the »PRO LOOP NX15 LOS« where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product.

When not used, store the device in a dry place, protected from dust.

Warranty

The »PRO LOOP NX7 LOS« is a very reliable product. Should a malfunction occur despite the unit being set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge.

It is recommended that you send in the product in its original packaging, so keep the packaging for the duration of the warranty period.

The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the product seal). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Disposal of used electric and electronic units (applicable in the countries of the European Union and other European countries with a separate collection system).

The symbol on the product or the packaging indicates that this product is not to be handled as ordinary household waste but has to be returned to a collecting point for the recycling of electric and electronic units.

You protect the environment and health of your fellow men by the correct disposal of this products. Environment and health are endangered by a faulty disposal.

Material recycling helps to reduce the consumption of raw material. You will receive further information on the recycling of this product from your local community, your communal disposal company or your local dealer.

Specifications

Height / Width / Depth:	43 mm x 430 mm x 290 mm (19" 1HU)
Weight:	3,8 kg
Power supply:	100 - 265 V AC 50 / 60 Hz
Passive protection:	Fuse 5x20 T5A
Power consumption in idle mode	17 W, at 230V AC with 2 connected 1 Ω -loops at room temperature after 30 minutes stabilizing time
Maximum power consumption:	2 x 224 W
Maximum input power:	< = 944 VA
Cooling system:	Fanless
Automatic Gain Control:	Speech-optimised, dynamic range: > 40 dB
Metal Loss Correction (MLC):	0 - 4 dB / octave
Operational range:	0°C - 45°C, < 2000 m above sea level

Loop output:

Max. current:	2x <7 A RMS
Max. tension:	2x 34 V peak
Loop resistance DC:	0,5 - 3,0 Ω
Connectors:	2x NL4 Speaker Twist, assigned 1+/2+ and 1-/2-
Frequency range:	80-6000 Hz (+/- 1,5 dB)
Distortion:	< 1% at rated power output, 1 kHz

Inputs:

Input 1	5-20 mV / 2 k Ω / 48V optional (MIC) 180 mV - 4 V / 10 k Ω (Line) Switchable between Mic und Line level, Combo-socket
Input 2	5-20 mV / 2 k Ω / 48V optional (MIC) 180 mV - 4 V / 10 k Ω (Line) Switchable between Mic und Line level, 3-point. Euroblock plug
Input 3	100V priority input, 2-point Euroblock plug
Link in	6,3mm jack socket, balanced

Outputs:

Headphones	3,5 mm stereo jack socket
Link connectors	6,35 mm jack socket, balanced
Fault status connector	3-point Euroblock plug
Loop connectors	NL4 Speaker Twist, assigned 1+ parallel to 2+ 1- parallel to 2-

This device complies with the following EC directives:

- 2017 / 2102 / EC RoHS-directive
- 2012 / 19 / EC WEEE-directive
- 2014 / 35 / EC Low voltage directive
- 2014 / 30 / EC Electromagnetic Compatibility

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at **www.humantechnik.com**.

Technical specifications subject to change without prior notice.

Humantechnik Service-Partner

D

Germany

Humantechnik GmbH

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Tel.: +49 (0) 76 21/ 9 56 89-0
Fax: +49 (0) 76 21/ 9 56 89-70
E-mail: info@humantechnik.com

CH

Switzerland

Humantechnik GHL AG

Rastatterstrasse 9
CH-4057 Basel

Tel.: +41 (0) 61/ 6 93 22 60
Fax: +41 (0) 61/ 6 93 22 61
E-mail: info@humantechnik.com

F/B

France
Belgium

SMS

Audio Electronique Sàrl

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Tel.: +33 (0) 3 89/ 44 14 00
Fax: +33 (0) 3 89/ 44 62 13
E-mail: sms@audiofr.com

NL

Netherlands

Hoorexpert BV

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Tel.: +31 (0) 3 45/ 63 23 93
Fax: +31 (0) 3 45/ 63 29 19
E-mail: info@hoorexpert.nl

GB

Great Britain

Sarabec Ltd

15 High Force Road
GB-Middlesbrough TS2 1RH

Tel.: +44 (0) 16 42/ 24 77 89
Fax: +44 (0) 16 42/ 23 08 27
E-mail: enquiries@sarabec.co.uk

**For other service-partners
in Europe please contact:**

Humantechnik Germany

Tel.: +49 (0) 76 21/ 9 56 89-0
Fax: +49 (0) 76 21/ 9 56 89-70
Internet: www.humantechnik.com
E-mail: info@humantechnik.com

HUMANTECHNIK